



Snowden-Pencer™ laparoscopic scissors



Snowden-Pencer™ Switch-Blade™ reposable scissors

Assemble in three simple steps:
thread, spread, click.



Thread
tip to shaft



Spread
the handle open



Click
to engage handle



Click-fit tip attachment for easy
assembly



Extended insulation helps prevent
the exposure of any conductive
metal surfaces

Cost-effective solution with a new blade each surgery

Switch-Blade reposable scissors with reusable handles and disposable scissor tips help provide superior cutting performance. The click-fit scissor tips make them easier to assemble by providing both audible and tactile feedback to confirm proper assembly.

Based on the next-generation Snowden-Pencer ring handle design, the Switch-Blade handle sets a new standard in comfort for laparoscopic instruments and accommodates both traditional precision-style grip or palm grip applications.



Ordering information

Catalog no.	Description	Qty
89-5103	Disposable Curved Metz scissor tips	10/box
89-5303	Disposable Curved Mini-Metz scissor tips	10/box
89-5115	Disposable Straight scissor tips	10/box
SP90-1053	Reusable handle and shaft, 36cm, insulated	1 each
SP90-1253	Reusable handle and shaft, 45cm, insulated	1 each

Snowden-Pencer™ reposable take-apart scissors

Easy assembly and disassembly

To assemble:

- Turn the rotation knob to align it with the electrosurgical post
- Press and hold the assembly/release buttons
- Slide the insert into the rotation knob opening until the nose key is seated in the opening
- Remove your fingers from the buttons, and insert the shaft until you hear an audible snap
- Actuate the handle to complete assembly

To disassemble:

- Press and hold the assembly/release buttons
- Pull the insert out from the handle.

Assembly/Release buttons

Nose key



Perform each surgery with a new blade and insulated shaft

New Snowden-Pencer reposable take-apart scissors provide fresh, sharp scissors and a new insulated shaft for each surgery. These features help optimize cutting performance, support sterility and advance product safety. Plus, they are easy to take apart and assemble.

Snowden-Pencer reposable take-apart scissors help reduce the need to clean and inspect the shaft and insulation. You can be confident in a fresh blade and a sterile, insulated shaft every time.



A clean cut

Our engineers specifically designed these scissors to withstand demanding surgical conditions. The scissors blades have gone through cutting tests on various surgical materials. We also tested the blades in mixed cold and hot cutting to simulate real surgical expectations. The scissors tips maintained their quality performance while withstanding aggressive surgery of more than 50 total cuts.

Ordering information

Catalog no.	Description		
SP95	Reposable scissors non-ratcheted handle		
Catalog no.	Scissors blade description	Shaft width	Shaft length
SP8301	Curved Metz	5 mm	36 cm
SP8401	Curved Metz	5 mm	45 cm
SP8302	Mini-Metz	5 mm	36 cm
SP8402	Mini-Metz	5 mm	45 cm
SP8303	Hook	5 mm	36 cm
SP8403	Hook	5 mm	45 cm

To learn more about Snowden-Pencer reusable laparoscopic scissors, contact your local sales representative, send an email to GMB-VM-EUROPE@carefusion.com or visit carefusion.com/snowdenpencer.

carefusion.com



CareFusion
1500 Waukegan Road
McGaw Park, IL 60085 USA



CareFusion France 309 S.A.S.
8 bis rue de la Renaissance
F-44110 CHATEAUBRIANT

©2013 CareFusion Corporation or one of its subsidiaries. All rights reserved. Snowden-Pencer and Switch-Blade are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. Medical devices class IIb according to Medical Devices Directive 93/42/EC - CE marked by TÜV Product Service GmbH (0123). LIT13006en-ED01-Mar2013



CareFusion