

Certus
Health™

Ports
Surgical Retrieval Bags
Laparoscopic Instruments
Laparoscopic Consumables

Ports

Optical Dilating

5mm



5mm Port with Dilating Trocar
BLT305
Box 20

11mm / 12mm



11mm / 12mm Port with Optical Dilating Trocar
BLT311
BLT312
Box 20



Optical Dilating Trocars
103mm Cannulas

Autoshielded

5mm



5mm Port with Autoshielded Trocar
BTR105
Box 20

10mm / 12mm



10mm / 12mm Port with Autoshielded Trocar
BTR110
BTR112
Box 20



Autoshielded Trocar
100mm Cannula

Hassan

10mm / 12mm



10mm / 12mm Hassan Trocar
BDS-TR03A
BDS-TR03B
Box 10

Retractor



A-Lap Curved Retractor Set
5mm - 6mm
MAS0210-01
Box 3

A-Lap™ is the first 3-Dimensional "hand-like" laparoscopic retractor. Features include a soft and compliant mesh controlled by an intuitive handle design, allowing the surgeon

to manipulate and lock-in the degree to which the device is opened. This can be used for General Surgery, Gynaecology and Urology procedures.

Retrieval Bags



10mm

Endo Retrieval Pouch
80mm - 200ml
BDS-D80
Box 10

10mm



Endo Retrieval Pouch
100mm - 350ml
BDS-D100
Box 10

10mm



Endo Retrieval Pouch
130mm - 700ml
BDS-D130
Box 10

Laparoscopic Instruments

5mm



Metz Scissor
BDS0010
Box 10

5mm



Babcock
BDS7500
Box 10

5mm



Endoclinch
BDS03-018
Box 10

5mm



Allis
BDS03-019
Box 10

5mm



Maryland
BDS02-002
Box 10

5mm



Short Johanne
BDS03-017
Box 10

5mm




Long Johanne
BDS03-021
Box 10

5mm



Extra Long Johanne for Bariatrics - 440mm Shaft
BDS-03-021-440
Box 10

Endoclinch
BDS03-018



Bipolar Marylands
BDS-0BIM/5

5mm



Bipolar Atraumatic Forceps
BDS0008
Box 5

5mm




Bipolar Marylands
BDS-0BIM/5
Box 5



Bipolar Cable
BDS0009
Box 25



Monopolar Cable
BDS1000
Box 25



Metz Scissor
BDS0010

Short Johanne
BDS03-017

Laparoscopic Consumables



Uterine Manipulator
BDS-UM
Box 20



Scope Warmer
BDS5000
Box 25

5mm



Suction Irrigation
BDS8010
Box 10

5mm



Suction Irrigation Double Spike
BDS8000
Box 12

5mm



Duel Position Suction Irrigation Double Spike with Laser Channel
BDS8020
Box 10



Disposable Knot Pusher
BDS-TP01
Box 10



Laparoscopic Smoke Filter
BDS2000
Box 25




Insufflation Tubing with Filter
BDS3000
Box 25

120mm / 150mm



Verres Needle 120mm
BO4-01-1005
Box 25
Verres Needle 150mm
BO4-01-1006
Box 25



Suture Grasper
BDS-SN01-C
Box 10



Ranfac Suture Grasper
RSG-14
Box 10



Anti-Fog Solution Bottle & Sponge
BDSAF01
Box 20



Camera Cover
BDS9010
Box 25



Adjustable Camera Drape
BDS9020
Box 25

Distributor

Contact

Certus Health™
Unit 7 Kingsford Rural Business Centre
Kentsbeare
Cullompton
Devon
EX15 2AU

T: 01884 266 377
F: 01884 266 388

www.certus-health.com