



Just when you think you've
seen everything.



GENIPalm™



BY THE SOCIETY OF
LAPAROENDOSCOPIC SURGEONS

Bladeless Tip Trocar and Cannula System

New Ideas for Modern Health Care

Access Ports: Palm Bladeless Tip Trocar and Cannula System



Ergonomic wing grip fits comfortably into palm while its low profile provides comfort and stability



Industry-first handle design, made from a single polycarbonate mold providing sturdiness and reliability with trocar placement



A distinctively engineered bladeless tip reduces insertion force to enter the abdominal cavity



Lean design reduces cost to produce product, savings passed on to the customer and a perfect option for budget-conscious healthcare systems



The uniquely-shaped bladeless tip separates and dissects tissue rather than cutting it, greatly reducing wound defect size

			One Cannula	Two Cannula
Diameter	Length	Qty/Box	Threaded	Threaded
5mm	70mm	5	100-005-051	100-005-052
5mm	100mm	5	100-005-151	100-005-152
10mm	100mm	5	100-010-151	---
12mm	100mm	5	100-012-151	---

Configurations with Smooth Cannula available for order in minimum quantities