



Sleek Vision all around



Natura Optical Bladeless Trocar

Access Ports: Optical Bladeless Trocar and Cannula System

The distinctively engineered visual tip greatly reduces wound defect size as well as insertion force into the abdominal cavity.

✓ Dual seal valve technology minimizes loss of insufflation pressure, maintaining pneumoperitoneum

✓ Injuries due to blind entry and a sharp blade are virtually eliminated

✓ Ergonomic wing handle provides superior comfort and stability

✓ Cannula incorporates a beveled and angled distal tip



High-visualization bladeless tip allows for visualization during insertion

Diameter	Length	Qty/Box	One Cannula	Two Cannula
			Threaded	Threaded
5mm	70mm	5	210-005-051	210-005-052
5mm	100mm	5	210-005-151	210-005-152
5mm	150mm	5	210-005-251	210-005-252
10mm	100mm	5	210-010-151	210-010-152
10mm	150mm	5	210-010-251	---
12mm	100mm	5	210-012-151	---
12mm	150mm	5	210-012-251	---
15mm	100mm	5	210-015-151	---
15mm	150mm	5	210-015-251	---

Every 15mm Cannula comes with a 5mm Vortex Reducer Valve