



Simple never looked so good.




Pyramidal Tip Trocar

New Ideas for Modern Health Care

Access Ports: Pyramidal Tip Trocar


The Geniport™ is designed to offer unrivaled versatility, efficiency, easy of use, and optimal cost.

Ergonomic handle design allows for a firm, controlled grip, while minimizing hand fatigue



Plastic-tipped Pyramidal Trocar
single-piece design

Plastic-tipped Pyramidal Trocar uses the same insertion force as a comparable metal-tipped trocar



5mm instrument diameter, in 70mm and 100 mm lengths



Available with smooth or threaded cannula



Available in several different configurations

			One Cannula	Two Cannula
Diameter	Length	Qty/Box	Threaded	Threaded
5mm	70mm	5	100-005-001	100-005-002
5mm	100mm	5	100-005-101	---

Configurations with Smooth Cannula available for order in minimum quantities