



Protection: on our mind
and built into our products.



Shielded Tip Trocar

Access Ports: Shielded Tip Trocar

The GeniShield™ is designed to offer unrivaled versatility, efficiency, ease of use, and optimal cost



The GeniShield™ is designed to be versatile, operationally flexible, and clinically adaptive to all surgical needs

Ergonomic handle design allows for a firm, controlled grip, while minimizing hand fatigue

Shielded Tip Trocar with retractable shield



One Cannula

Two Cannula

| Diameter | Length | Qty/Box | Threaded | Threaded |
|----------|--------|---------|-------------|-------------|
| 5mm | 70mm | 5 | 100-005-091 | 100-005-092 |
| 5mm | 100mm | 5 | 100-005-191 | 100-005-192 |
| 10mm | 100mm | 5 | 100-010-191 | 100-010-192 |
| 12mm | 100mm | 5 | 100-012-191 | --- |

Configurations with Smooth Cannula available for order in minimum quantities

6869 Stapoint Court #114 | Winter Park, Florida USA 32792 | Tel: 1.407.657.4851 Fax: 1.407.677.9773 | www.geniconendo.com

© 1998-2016 GENICON® All rights Reserved. US and foreign patents issued or pending.
Visit www.geniconendo.com for a list of trademarks and patents

Rev. D