



# CINCH<sup>®</sup>

## Tube Secure Device

**NEW! The CINCH<sup>®</sup> Family Now Includes Two Sizes!**



**Medium CINCH<sup>®</sup>**  
4F to 17F



**Large CINCH<sup>®</sup>**  
9F to 30F

### Patient Comfort First

The low profile and flexible design allows for comfort and mobility.

### Safe + Secure

Anchors medical tubes while relieving stress and strain at insertion points.

### Hypoallergenic Adhesive

Provides gentle yet secure attachment.

### Easy to Use

Simply clean the skin, apply the CINCH<sup>®</sup>, and secure the tube!



# CINCH®


From 4F to 30F, the CINCH® Has You Covered!

## Versatile

Medium CINCH® secures tubes 4F to 17F,  
Large CINCH® secures tubes 9F to 30F.

## Comfortable

Soft, low-profile, flexible.



## Breathable Backing


In two convenient sizes,  
Medium (3 cm x 4 cm) and  
Large (4 cm x 6 cm).

## Hypoallergenic

Gentle adhesive for sensitive skin.

## Reusable


Elastomer strap can be  
repeatedly opened and closed.




The CINCH® comes fully assembled with an  
easy-peel split backing.

**CINCH®**

## Medium CINCH®



## Large CINCH®



Description	Tube Securement Range	Box Qty	Order Number
Medium CINCH®	4F to 17F	10	CINCH417M-10
		50	CINCH417M-50
Large CINCH®	9F to 30F	10	CINCH930L-10
		50	CINCH930L-50

For a complete list of products, visit us at <https://www.kebomed.co.uk/products/amt-enteral-devices/>

01884 266377 | [enteral.sales@kebomed.co.uk](mailto:enteral.sales@kebomed.co.uk)

# KEBOMED

Kebomed UK  
Unit 7, Kingsford Rural Business Centre  
Kentisbeare  
Cullompton EX15 2AU

© 2021 Applied Medical Technology, Inc. C-4733-A

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>  
CAUTION: If redness or irritation occurs, stop use and consult a healthcare professional. Do not over tighten.  
Images not to scale. This document is not a substitute for the DFU (Directions for Use).