


Comprehensive I.P. Kits

Initial Placement Gastrostomy Sets

The intuitive design of AMT I.P. Kits ensures all of the essential tools for initial placement of a percutaneous gastrostomy feeding device are packaged in one easy-to-use kit.



Quality, Right from the Start

- 1 Introducer Needles**
Two lengths (short, long) to accommodate a broad range of patients.
- 2 Introducer / Stiffener**
An AMT exclusive, size-specific, guidewire compatible.
- 3 Slip Tip Syringe**
For use with Balloon Stoma Measuring Device.
- 4 BSMD**
Guidewire compatible Balloon Stoma Measuring Device.
- 5 Safety Scalpel**
- 6 Gauze**
- 7 Guidewire**
With J-tip
- 8 Multiple Size Dilators**

Designed to Minimize Trauma Associated with Stoma Tract Development

- Tapered Dilators Help Ensure Proper Dilation
- Requires Only One Endoscopic Pass
- Compatible with Laparoscopic, Endoscopic, and Fluoroscopic Placement
- For Use with 10F - 18F Devices



Innovating. Educating. Changing Lives.™

Comprehensive I.P. Kits

Intuitive Design. AMT Quality.

AMT offers Comprehensive and Basic I.P. Kits that are designed to place 10F - 18F gastrostomy devices. Comprehensive kits include dilators in a range of French sizes, a guidewire, two introducer needles (short, long), a balloon stoma measuring device, an introducer/stiffener, a scalpel, a small syringe, and gauze. Basic kits include dilators in a range of French sizes, a guidewire, two introducer needles (short, long), and an introducer/stiffener.

Comprehensive I.P. Kits[‡]

ORDER #	KIT CONTENTS
IP-10-S	(3) Dilators (7F, 12F, 14F-S)
	(1) Guidewire, 70 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(1) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
(4) 4" x 4" Gauze	
ORDER #	KIT CONTENTS
IP-1218	(5) Dilators (7F, 12F, 16F, 18F, 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
(4) 4" x 4" Gauze	
ORDER #	KIT CONTENTS
IP-1218-S	(5) Dilators (7F, 12F, 16F, 18F-S , 22F)
	(1) Guidewire, 90 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible BSMD
	(2) Guidewire Compatible Introducer/Stiffener
	(1) Safety Scalpel, #11
	(1) 5 cc Syringe
(4) 4" x 4" Gauze	

I.P. Kits are packaged 1/box (sterile). | [‡]NOTE: Kits do not contain T-Fasteners or U-Stitches.
^{-S} I.P. Kits include one sheathed dilator, noted in bold above.







Basic I.P. Kits^{*}

ORDER #	KIT CONTENTS
IP-DIL-10	(3) Dilators (7F, 12F, 14F)
	(1) Guidewire, 45 cm (0.018")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener
ORDER #	KIT CONTENTS
IP-DIL	(4) Dilators (7F, 12F, 16F, 20F)
	(1) Guidewire, 45 cm (0.038")
	(2) Introducer Needles, 18 Gauge (short, long)
	(1) Guidewire Compatible Introducer/Stiffener

I.P. Kits are packaged 1/box (sterile). | * NOTE: Kits do not contain T-Fasteners or U-Stitches.

I.P. Kit Accessories

Available Separately

ORDER #	DESCRIPTION
1804	Balloon Stoma Measuring Device Guidewire Compatible (up to 0.038") 
1704	Balloon Stoma Measuring Device Non-Guidewire Compatible 
4-6010	Introducer/Stiffener Guidewire Compatible (up to 0.018") Device Compatibility: 10F and larger Low Profile Devices 
4-6000	Introducer/Stiffener - 4.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices 
4-6200	Introducer/Stiffener - 9.6" Length Guidewire Compatible (up to 0.038") Device Compatibility: 12F and larger Low Profile Devices 
4-5000	Introducer/Stiffener Non-Guidewire Compatible Device Compatibility: 12F and larger MINI™ Button Devices 

I.P. Kit accessories are individually pouched and sold 10/box (sterile).

For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



Comprehensive I.P. Kits



Applied Medical Technology, Inc.
 8006 Katherine Boulevard
 Brecksville, OH 44141
www.AppliedMedical.net
 An ISO 13485 Registered Company
 Cleared For Initial Placement
 Made In USA

