



Urology

Uromedex
Endourology range

MEDICAL DEVICE SOLUTIONS

UROMED

Your Competent Partner for Urological Medical Devices

over more than 50 years, UROMED has developed into one of the leading providers of urological medical devices. UROMED offers a comprehensive range of products for diagnosis and therapy of urological diseases at Oststeinbek, close to Hamburg. At UROMED, a team of around 190 employees ensures that the family business, founded in 1967, continuously develops furthermore and keeps the outstanding reputation as a specialist in urology.

Based on long-term experience and know-how, UROMED offers its customers individual solutions to meet the requirements of users and consumers. They develop appropriate solutions using proven as well innovative technologies that are sustainable and economically viable.



TABLE OF CONTENTS

Transurethral Catheters®	4
Post Operative Catheters	5
Suprapubic Urinary Catheters	6
Ureteral Catheter	7
Heparius 12 month indwelling ureteral stent	8
Hydropur 6 month indwelling ureteral stent	9
Catheterisation accessories	10
UROMED STONIZER® advanced Stone Retrieval Basket	11

Transurethral Catheters®

biocompatible pure silicone for the required safety, shape and consistency for comfort and long-term application



Whistle Tip



Tiemann Tip



Silucent Tieman Tip


Transparent and pure silicone



Short Neleton Tip



Neleton Tip



Silucent Nesleton Tip

Transparent and pure silicone

Order Information

DESCRIPTION	ORDER #	BALLOON SIZE	BOX/QTY
Whistle Tip Catheter			
Short Open Whistle tip, 1 eye	1091	3-10 ml	10
Open Whistle tip, 1 eye	1090	30 ml	10
Tiemann Tip Catheter			
Tiemann tip, 1 eye	1240	3-10 ml	10
Olive Tiemann tip, 1 eye, soft	1245	3-10	10
Silucent Tiemann Tip Catheter			
Olive Tiemann tip, 1 eye, Transparent	1250	3-10	10
Short Neleton Tip Catheter			
short, conical tip, 2 eyes	1095	3-10 ml	10
Short, conical tip,	1100	3-10 ml	10
Short, conical tip, central open application over guidewire	1111	5-10 ml	10
Short, conical tip, integrated flat balloon	1110	5-10 ml	10
Pediatric Catheter, short conical tip, plastic guidewire	1340	1-3 ml	10
Female catheter, short conical tip	1370	3-10 ml	10
Neleton Tip Catheter			
Conical tip, 2 eyes	1101	30 ml	10
Silucent Neleton Tip Catheter			
Short, conical tip, transparent	1104	3-10 ml	10
Conical tip, transparent, 2 eyes	1105	30 ml	10

Post Operative Catheters

A reliable catheter ensures the postoperative success of TUR and other surgeries. Our efficient silicone catheters are the result of a continuous demand from surgeons for improvement of OP catheters.



Continuous Irrigation



Continuous Irrigation
with central hole



DBK Couvelaire



DBK Delionette



DBK Dufour



Neobladder by Hanke



Neobladder DBK



TBK-Delionette



TBK-Dufour

Order Information

DESCRIPTION	ORDER #	BALLOON SIZE	BOX/QTY
Conituous Irrigation Catheter			
Cylindrical tip, 3-way, 2 eyes, length 43 cm	2080	30-50 ml	10
Continuous Irrigation with Central Hole Catheter			
Cylindrical tip, 3-way, 2 eyes, length 43 cm	2085	40-50	10
DBK Couvelair			
Short, nonkinking Couvelaire tip 3-way, 2 eyes, l length 43 cm, solid	2130	40-50	10
Short, nonkinking Couvelaire tip 3-way, 2 eyes, l length 43 cm, soft	2132	40-50	10
DBK Delionette			
Long, nonkinking Delinotte tip 3-way, 2 eyes, length 43 cm	2240	30-70 ml	10
DBK Dufour			
long, nonkinking Dufour tip 3-way, 2 eyes, length 43 cm	2340	40-70 ml	10
Neobladder by Hanke			
5 cm long, central open tip 3-way, 4 large eyes, length 45 cm	3220	30 ml	10
Neobladder DBK			
5 cm long, central open tip 4 large eyes, length 45 cm	3221	30 ml	10
TBK Delinotte			
Delinotte tip 2-way, 2 eyes, length 43 cm	1770	40-60 ml	10
TBK Dufour			
Dufour tip 2-way, 2 eyes, length 43 cm	1870	40-60 ml	10

Suprapubic Urinary Catheters

Central to the Suprapubic range is the incorporation of a balloon catheter which makes practised fixation at the skin unnecessary. Further decisive steps like INTEGRAL-, SUPRASOFT- and DILATION CATHETERS enrich our program and continue to meet the demands for a safe and comfortable application.



**UROMED SUPRA-KATH
INTEGRAL**




**UROMED Suprapubic
Balloon Catheter**



**UROMED »SILUCENT«
SUPRA-KATH-INTEGRAL**



**UROMED »SUPRAPUR«
puncture Set in
Seldinger-technique**



**UROMedCYST®
»INTEGRAL« Balloon
Catheter Exchange Set**




**UROMED »SUPRASOFT«
Double Balloon Catheter**

Order Information


DESCRIPTION	ORDER #	BALLOON SIZE	BOX/QTY
UROMED SUPRA-KATH INTEGRAL			
Short, atraumatic and central open tip with rounded edge integrated flat balloon 2 eyes, length 33 cm	2965	5 ml	10
Long, atraumatic and central open tip with rounded edge Integrated flat balloon 2 eyes, length 43 cm	2955	5 ml	10
UROMED Suprapubic Balloon Catheter			
Short, atraumatic and central open tip with rounded edge 2 eyes, length 33 cm	2966	2-5 ml	10
Long, atraumatic and central open tip with rounded edge 2 eyes, length 43 cm	2956	2-5 ml	10
UROMED »SILUCENT« SUPRA-KATH-INTEGRAL			
atraumatic and central open tip with rounded edge 2 eyes, length 43 cm transparent	2953	5-10 ml	10
UROMED »SUPRAPUR« puncture Set in Seldinger-technique			
<ul style="list-style-type: none"> dilation balloon catheter »INTEGRAL« cylindrical open tip, length 43 cm (REF 3045.12) Lunderquist guidewire with insertion guide 0,035"/70 cm puncture cannula 17,5 GA/1,3 mm, needle length 18,5 cm Fascial dilator Ch. 10/3,3 mm (REF 2650) Fascial dilator Ch. 14/4,7 mm (REF 2650) catheter plug 	3047	5 ml	10
UROMedCYST® »INTEGRAL« Balloon Catheter Exchange Set			
short, atraumatic tip with rounded edge x-ray contrast line	3021	5 ml	10
UROMED »SUPRASOFT« Double Balloon Catheter			
without tip, central open 1 eye, length 36 cm especially for tenesmen	2960	15-20 ml	10

Ureteral Catheter

Uromed open ureteral catheters offer instruments with an extremely high gliding ability and wide lumen for application in Seldinger-technique. Additionally they also permit the use of thicker guidewires



**UROMED PURgreen®
Ureter Catheter**



**UROMED PURgreen®
Ureter Catheter
by Chevassu**



**UROMED Ureteral
Catheter**




**UROMED URETERAL
CATHETER »DUOKATH**

Order Information

DESCRIPTION	ORDER #	Channel size	BOX/QTY
UROMED PURgreen® Ureter Catheter			
Atraumatic eye, length 70 cm Cylindrical tip high gliding ability	4070	Ch.3-8	10
Atraumatic eye, length 70 cm Tiemann tip high gliding ability	4080	Ch. 3-7	10
UROMED PURgreen® Ureter Catheter by Chevassu			
Catheter for the instillation of a contrast medium into the ureter Button-whistle-tip atraumatic eye, length 70 cm	4090	Shaft/Tip Ch. 3/Ch. 5 Ch. 4/Ch. 6 Ch. 4/Ch. 7 Ch. 5/Ch. 8	10
UROMED Ureteral Catheter			
central open Tiemann tip length 70 cm	4081	Ch.5	10
UROMED Ureter Catheter Duocath			
<ul style="list-style-type: none"> Catheter for retrograde probing of the ureter 2 in 1! The ureteral catheter with different tips for application on both sides and many advantages. central open Tiemann tip and conical tip length 70 cm 	4082	Ch.5	10

Heparius 12 month indwelling ureteral stent


Uromed open ureteral catheters offer instruments with an extremely high gliding ability and wide lumen for application in Seldinger-technique. Additionally they also permit the use of thicker guidewires



**UROMED HEPARIUS®
Ureteral Stent**



**UROMED HEPARIUS®
Ureteral Stent –
Fixing clamp**




**UROMED HEPARIUS®
Ureteral Stent
– Guidewire**

Order Information

DESCRIPTION	ORDER #	CHANNEL SIZE	BOX/QTY
UROMED HEPARIUS® Ureteral Stent			
• Ureteral stent, double pigtail	4145	Ch.4.8	10
• Central open tip,	4146	Ch.6.0	10
• Perforation in the ureteral part	4147	Ch.7.0	10
• Pusher length, Ch 4.8 = 70cm Ch. 6 + 7 = 45 cm			
• Ureteral stent, double pigtail	4165	Ch.4.8	10
• Central open tip,	4166	Ch.6.0	10
• No Perforation in the ureteral part	4167	Ch.7.0	10
• Pusher length, Ch 4.8 = 70cm Ch. 6 + 7 = 45 cm			
UROMED HEPARIUS® Ureteral Stent with fixing clamp			
• Ureteral stent, double pigtail	4175	Ch.4.8	10
• Central open tip,	4176	Ch.6	10
• Perforation in the ureteral part	4177	Ch.7	10
• Pusher length, Ch 4.8 = 70cm Ch. 6 + 7 = 45 cm			
UROMED HEPARIUS® Ureteral Stent with guidewire			
• Ureteral stent, double pigtail	4155	Ch.4.8	10
• Central open tip,	4156	Ch.6	10
• Perforation in the ureteral part	4157	Ch.7	10
• Pusher length, Ch 4.8 = 70cm Ch. 6 + 7 = 45 cm			
• PTFE coated guidewire			

Hydropur 6 month indwelling ureteral stent


HYDROPUR®« convinces by its high gliding ability, stability during the long indwelling time and easy handling. The stents are available in finely graduated sizes as well as manoeuvrable versions or with pusher an application thread.




**UROMED HYDROPUR®
Ureteral Stent**



**UROMED HYDROPUR®
Ureteral Stent
Fixing clamp set**



**UROMED HYDROPUR®
Ureteral Stent
Guidewire set**




**UROMED HYDROPUR®
Ureteral Stent External
drainage set**

Order Information

DESCRIPTION	ORDER #	CHANNEL SIZE	BOX/QTY
UROMED HEPARIUS® Ureteral Stent			
• central open tip	4455	Ch.4.8	10
• without perforation in the ureteral part	4456	Ch.6.0	10
• For stents Ch. 4,8: guidewire 0,028"	4457	Ch.7.0	10
• For stents Ch. 6-8: guidewire 0,035"	4458	Ch.8.0	10
• central open tip	4445	Ch.4.8	10
• with perforation in the ureteral part	4446	Ch.6.0	10
• For stents Ch. 4,8: guidewire 0,028"	4447	Ch.7.0	10
• For stents Ch. 6-8: guidewire 0,035"	4448	Ch.8.0	10
UROMED HEPARIUS® Ureteral Stent with fixing clamp			
• The ureteral stent in Ch. 4,8 is for ureterorenoscopic application.	4435	Ch.4.8	10
	4436	Ch.6	10
• central open tip	4437	Ch.7	10
• with perforation in the ureteral part	4438	Ch.8	10
	4439	Ch.9	10
UROMED HEPARIUS® Ureteral Stent with fixing clamp and guidewire			
• The ureteral stent in Ch. 4,8 is for ureterorenoscopic application.	4405	Ch.4.8	10
	4406	Ch.6	10
• cylindrical tip	4407	Ch.7	10
• with perforation in the ureteral part	4408	Ch.8	10
• PTFE Guidewire	4409	Ch.9	10
UROMED HYDROPUR® Ureteral Stent External drainage set			
• cylindrical tip	4465	Ch.4.8	10
• without perforation in the ureteral part	4466	Ch.6	10
	4467	Ch.7	10
• aromatic polyurethane	4478	Ch.8	10
• With Fixation clamp			
• central open tip	4475	Ch.4.8	10
• without perforation in the ureteral part	4476	Ch.6	10
	4477	Ch.7	10
• aromatic polyurethane	4478	Ch.8	10
	4479	Ch.9	10

Catheterisation accessories

With the UROMED catheter valve the patient can counteract methodically and promote the pre-servation of natural bladder function.



**UROMED basic
Catheter Valve**



**UROMED Compact
Catheter Valve**




**UROMED Universal
Catheter Valve**

Order Information


DESCRIPTION	ORDER #	BOX/QTY
UROMED Basic Catheter Valve		
<ul style="list-style-type: none"> permanent opening possible by lever repeated insertion and detachment of a catheter plug is not required 	1506	5
UROMED Compact Catheter Valve		
<ul style="list-style-type: none"> should be used if the connection with an urine bag is not required repeated insertion and detachment of a catheter plug is not required one-handed usability of the push-button 	1504	5
UROMED Universal Catheter Valve		
<ul style="list-style-type: none"> UROMED catheter valve 	1500	5
<ul style="list-style-type: none"> UROMED catheter valve with catheter valve adapter 	1501	5
<ul style="list-style-type: none"> UROMED catheter valve adapter with safety bayonet lock for connection of a catheter valve with a UROMED urine bag with step connector 	1505	5

UROMED STONIZER® advanced Stone Retrieval Basket


Retrieval basket for the urological removal of stones during ureterorenoscopy



Tipless, 4-wire
basket, straight,
3-finger-handle




UROMED STONIZER®
ROTAR360



With tip, 4-wire basket,
straight/helical,
3-finger handle



With tip, 4-wire
basket, straight/
helical, stick handle



Open 3-wire grasper,
stick handle

Order Information

DESCRIPTION	ORDER #	Ch./Fr	Ø Shaft mm	Ø Basket mm
Tipless, 4-wire basket, straight, 3-finger-handle				
Straight basket	78259016	2.5	0.8	16
Straight basket	78309016	3.0	1.0	16
Helical Basket	79309016	3.0	1.0	16
Helical basket	79409016	4.0	1.3	16
UROMED STONIZER® ROTAR360				
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71151209	1.5	0.5	9
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71151212	1.5	0.5	12
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71181209	1.8	0.6	9
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71181212	1.8	0.6	12
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71181215	1.8	0.6	15
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71221212	2.2	0.7	12
UROMED STONIZER® ROTAR360 Stone Retrieval Basket	71221215	2.2	0.7	15
With tip, 4-wire basket, straight/helical, 3-finger handle				
Straight basket	78259016	2.5	0.8	16
Straight basket	78309016	3.0	1.0	16
Helical Basket	79309016	3.0	1.0	16
Helical basket	79409016	4.0	1.3	16
With tip, 4-wire basket, straight/helical, stick handle				
Straight basket	74259015	2.5	0.8	15
Straight basket	74309015	3.0	1.0	15
Helical Basket	75309015	3.0	1.0	15
Helical basket	75409015	4.0	1.3	15
Open 3-wire grasper, stick handle				
UROMED STONIZER® TWONE® Stone Extractor	72151208	1.5	0.5	8
UROMED STONIZER® TWONE® Stone Extractor	72151211	1.5	0.5	11
UROMED STONIZER® TWONE® Stone Extractor	72201208	2.0	0.6	8
UROMED STONIZER® TWONE® Stone Extractor	72201211	2.0	0.6	11