

MiniONE Button

Capsule Non-Balloon Devices



Exclusive Capsule Technology

Material facilitates easy and comfortable insertion.

Glow Green™ Technology

Easier nighttime feeds.

Improve Comfort

Combines the advantage of a low-profile device with the longer life of a non-balloon device.

Non-Balloon Advantage

Ideal for patients who don't tolerate balloon buttons.

More Durable

Internal bolster is smaller than most devices, durable and less likely to block the pylorus.¹



Pull for Capsule Deployment



MiniONE Button

Capsule Non-Balloon Devices

MiniONE®

Order Information

	LENGTH (cm)	ENFit® ORDER 14Fr	ENFit® ORDER 18Fr	ENFit® ORDER 20Fr	ENFit® ORDER 24Fr
Capsule Non-Balloon Button	1	M1-2-1410-I			
	1.2	M1-2-1412-I	M1-2-1812-I		
	1.5	M1-2-1415-I	M1-2-1815-I		
	1.7	M1-2-1417-I	M1-2-1817-I	M1-2-2017-I	M1-2-2417-I
	2	M1-2-1420-I	M1-2-1820-I		
	2.4		M1-2-1824-I	M1-2-2024-I	M1-2-2424-I
	2.5	M1-2-1425-I			
	3	M1-2-1430-I	M1-2-1830-I		
	3.4		M1-2-1834-I	M1-2-2034-I	M1-2-2434-I
	4.4		M1-2-1844-I	M1-2-2044-I	M1-2-2444-I

Feed Sets Included in Kit



MiniONE® Button comes packaged as 1/box (sterile). Kit contents include device (1), bolus extension set (1), y-port extension set (1), 10cc syringe (1), 35cc bolus syringe (1), Placement T-Handle/Removal T-Handle†, snap arm assembly, and gauze.

Spare Removal Tools for Obturated Buttons

FRENCH SIZE	DESCRIPTION	ORDER #	BOX/QTY
14F	1.0, 1.2 cm	R-S1012	1
14F	1.5, 1.7 cm	R-S1517	1
14F	2.0 cm	R-S2023	1
14F	2.5 cm	R-S2427	1
14F	3.0 cm	R-S3033	1
18F	1.2 cm	R-L1012	1
18F, 20F, 24F	1.5, 1.7 cm	R-L1517	1
18F	2.0 cm	R-L2023	1
18F	2.4 cm	R-L2427	1
18F	3.0 cm	R-L3033	1
18F, 20F, 24F	3.4 cm	R-L3437	1
18F, 20F, 24F	4.4 cm	R-L4044	1

MT Obturator for Distending button during insertion and removal.



Capsule Non-Balloon Kit: Removal T-Handle

