

# Traditional G-JET<sup>®</sup>

Adjustable Length Gastric-Jejunal Enteral Feeding Tube



## Reduce Clogging

Unique internal structure transitions from a triple lumen to a single lumen, designed to reduce clogging.<sup>1</sup>

## Anti-Kink Technology\*

Spring reinforcement stabilizes the jejunal lumen to help reduce kinking and clogging.

## AMT Exclusive Balloon

Apple-shaped, expands horizontally creating a superior gastric seal.<sup>2</sup>

## Improve Comfort

The device is made from medical-grade silicone to reduce skin irritation and improve patient comfort.


The G-JET<sup>®</sup> You Trust, Now Offered in Traditional Length



KEBOMED

# Same French Size, Larger Lumens

DESIGNED TO REDUCE CLOGGING



## GASTRIC Lumen Comparison<sup>1</sup>

AMT 14F is **Exclusive**, Not Offered by Competitor

AMT 16F is **41% Larger** than Competitor 16F

AMT 18F is **34% Larger** than Competitor 18F

## JEJUNAL Lumen Comparison<sup>1</sup>

AMT 14F is **56% Larger** than Competitor 16F  
**22% Larger** than Competitor 18F

AMT 16F is **131% Larger** than Competitor 16F  
**81% Larger** than Competitor 18F

AMT 18F is **172% Larger** than Competitor 18F  
**57% Larger** than Competitor 22F

## Traditional Length G-JET® | Order Information\*

JEJUNAL	GASTRIC 14F Without Spring	GASTRIC 16F* With Spring	GASTRIC 18F* With Spring	GASTRIC 22F Without Spring
LENGTH (cm)	ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #	ENFIT® ORDER #
15	TGJ-14-15-I	TGJ-16-15-I		
22	TGJ-14-22-I	TGJ-16-22-I	TGJ-18-22-I	
30	TGJ-14-30-I	TGJ-16-30-I	TGJ-18-30-I	TGJ-22-30-I
45	TGJ-14-45-I	TGJ-16-45-I	TGJ-18-45-I	TGJ-22-45-I

Bal. Fill Vol.	MIN.	REC.	MAX.	MIN.	REC.	MAX.	MIN.	REC.	MAX.	MIN.	REC.	MAX.
	3 ml	4 ml	5 ml	3 ml	4 ml	5 ml	7 ml	8 ml	9 ml	8 ml	10 ml	12 ml

Traditional G-JET® comes packaged as 1/box (sterile). | Kit contents include: device (1), 5cc/10cc balloon syringe (1), 35cc bolus syringe (1), and gauze.

\*14F + 22F G-JET® Devices are MR Safe – 16F + 18F G-JET® Devices are MR Conditional.

NOTE: G-JET® devices should never be cut to custom lengths.

For a complete list of products, visit us at <https://www.kebomed.co.uk/products/amt-enteral-devices/>

01884 266377 | [enteral.sales@kebomed.co.uk](mailto:enteral.sales@kebomed.co.uk)

# KEBOMED

**KeBomed UK**  
Unit 7, Kingsford Rural Business Centre  
Kentisbeare  
Cullompton EX15 2AU

© 2024 Applied Medical Technology, Inc. C-4700-H

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

1.Claim based on internal testing of 14F-22F devices. | 2.Claims based on internal testing of 16F, 18F, and 22F devices. | ENFIT® is a registered trademark of GEDSA.