



# Urology

Marflow AG  
Soodstrasse 57  
8134 Adliswil/Zurich  
Switzerland  
Tel. +41 44 709 01 01  
Fax +41 44 710 21 36  
[marflow@marflow.ch](mailto:marflow@marflow.ch)  
[www.marflow.ch](http://www.marflow.ch)

**MarFlow**<sup>®</sup>  
SWITZERLAND


Index

| Page | Article                                           |
|------|---------------------------------------------------|
| 4    | <b>Ureter Stents – Standard</b>                   |
| 4    | <b>Ureter Stents – Multiloop</b>                  |
| 6    | <b>Ureter Stents – Longlife 180</b>               |
| 6    | <b>Ureter Stents – Longlife 360</b>               |
| 7    | <b>Ureter Stents – Groove</b>                     |
| 7    | <b>Ureter Stents – Hydrophilic coated</b>         |
| 8    | <b>Ureter Stents – Tumor</b>                      |
| 9    | <b>Ureter Stents – Tumor 100% Silicone</b>        |
| 10   | <b>Ureter Stents – Paediatric</b>                 |
| 10   | <b>Ureter Stents – Renal Transplant</b>           |
| 10   | <b>Ureter Stents – Intra Operative</b>            |
| 12   | <b>Ureter Stents – Mono-J</b>                     |
| 12   | <b>Blue Stents via Ureterorenoscope</b>           |
| 13   | <b>Steerable Pusher</b>                           |
| 14   | <b>Grasping Forceps flexible</b>                  |
| 15   | <b>Grasping Forceps metallic-stiff</b>            |
| 16   | <b>Stainless Steel Guidewires</b>                 |
| 16   | <b>Nitinol Super Guidewires</b>                   |
| 17   | <b>Cobra Nitinol Guidewires</b>                   |
| 18   | <b>Ureter Catheters</b>                           |
| 20   | <b>Ureter Access Sheath Standard</b>              |
| 21   | <b>Ureter Access Sheath Braided</b>               |
| 22   | <b>Kidney Stone Basket</b>                        |
| 23   | <b>MG-Paw Stone Basket</b>                        |
| 24   | <b>X-Circle Stone Basket</b>                      |
| 25   | <b>Kidney Stone Stop</b>                          |
| 26   | <b>Nephrostomy Drainage Set PIGTAIL</b>           |
| 28   | <b>Nephrostomy Drainage Set MALECOT</b>           |
| 30   | <b>Nephrostomy Drainage Set INTEGRAL</b>          |
| 32   | <b>Nephrostomy Drainage Set Locking Mechanism</b> |
| 33   | <b>Nephrostomy Drainage Set Direct Puncture</b>   |
| 34   | <b>Re-Entry MALECOT Drainage Catheter</b>         |
| 35   | <b>Initial Puncture Needles</b>                   |

| Page | Article                                     |
|------|---------------------------------------------|
| 36   | <b>Fascial Dilators</b>                     |
| 37   | <b>Ureter Dilator Set</b>                   |
| 38   | <b>Ureter Balloon-Dilator</b>               |
| 38   | <b>Transureteroscopic Balloon-Dilator</b>   |
| 38   | <b>Nephrostomy Track Balloon-Dilator</b>    |
| 39   | <b>Inflation Device</b>                     |
| 40   | <b>Nottingham One-Step Dilator</b>          |
| 40   | <b>Female and Meatal Dilator</b>            |
| 41   | <b>Urethra Dilator Set Curved</b>           |
| 42   | <b>Amplatz Renal Dilation Set</b>           |
| 43   | <b>Amplatz Renal Dilators + Sheath</b>      |
| 44   | <b>Mini PCNL – Dilation Set</b>             |
| 45   | <b>Semi-Automatic Biopsy-Gun</b>            |
| 46   | <b>Suprapubic Drainage CYSTOCATH</b>        |
| 48   | <b>Suprapubic Drainage INTEGRAL</b>         |
| 49   | <b>Suprapubic Drainage PIGTAIL</b>          |
| 50   | <b>Haematuria Catheter</b>                  |
| 52   | <b>TUR Electrodes</b>                       |
| 53   | <b>TUR Electrodes MONOPOLAR</b>             |
| 54   | <b>TUR Electrodes BIPOLAR</b>               |
| 55   | <b>Bladder Evaquator</b>                    |
| 56   | <b>Bladder Irrigation Set</b>               |
| 57   | <b>Urine Bag</b>                            |
| 58   | <b>2-Way Balloon Catheter</b>               |
| 60   | <b>Catheter Gel</b>                         |
| 60   | <b>Sterile Liquid with 10% Glycerine</b>    |
| 61   | <b>Hydrophilic Urinary Catheters + Kits</b> |
| 62   | <b>Intra-Uterine Insemination-Catheter</b>  |
| 63   | <b>Dialyse Catheter</b>                     |

Ureter Stents  
- Standard Double-J  
- Multiloop

Made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days



Ureter Stents

Set consists of stent, pusher, suture and clamp. Made of Polyurethane Tecoflex® USA

Standard Double-J

| One side open  | Both sides open | Size (Fr/CH) | Length (cm) | Guidewire |
|----------------|-----------------|--------------|-------------|-----------|
| <b>SS-4716</b> | <b>SOT-4716</b> | 4.7          | 16          | 0.032"    |
| <b>SS-4724</b> | <b>SOT-4724</b> | 4.7          | 24          | 0.032"    |
| <b>SS-4726</b> | <b>SOT-4726</b> | 4.7          | 26          | 0.032"    |
| <b>SS-4728</b> | <b>SOT-4728</b> | 4.7          | 28          | 0.032"    |
| <b>SS-4730</b> | <b>SOT-4730</b> | 4.7          | 30          | 0.032"    |
|                |                 |              |             |           |
| <b>SS-516</b>  | <b>SOT-516</b>  | 5            | 16          | 0.035"    |
| <b>SS-524</b>  | <b>SOT-524</b>  | 5            | 24          | 0.035"    |
| <b>SS-526</b>  | <b>SOT-526</b>  | 5            | 26          | 0.035"    |
| <b>SS-528</b>  | <b>SOT-528</b>  | 5            | 28          | 0.035"    |
| <b>SS-530</b>  | <b>SOT-530</b>  | 5            | 30          | 0.035"    |
|                |                 |              |             |           |
| <b>SS-616</b>  | <b>SOT-616</b>  | 6            | 16          | 0.035"    |
| <b>SS-624</b>  | <b>SOT-624</b>  | 6            | 24          | 0.035"    |
| <b>SS-626</b>  | <b>SOT-626</b>  | 6            | 26          | 0.035"    |
| <b>SS-628</b>  | <b>SOT-628</b>  | 6            | 28          | 0.035"    |
| <b>SS-630</b>  | <b>SOT-630</b>  | 6            | 30          | 0.035"    |
|                |                 |              |             |           |
| <b>SS-716</b>  | <b>SOT-716</b>  | 7            | 16          | 0.038"    |
| <b>SS-724</b>  | <b>SOT-724</b>  | 7            | 24          | 0.038"    |
| <b>SS-726</b>  | <b>SOT-726</b>  | 7            | 26          | 0.038"    |
| <b>SS-728</b>  | <b>SOT-728</b>  | 7            | 28          | 0.038"    |
| <b>SS-730</b>  | <b>SOT-730</b>  | 7            | 30          | 0.038"    |
|                |                 |              |             |           |
| <b>SS-826</b>  | <b>SOT-826</b>  | 8            | 26          | 0.038"    |
| <b>SS-828</b>  | <b>SOT-828</b>  | 8            | 28          | 0.038"    |
| <b>SS-830</b>  | <b>SOT-830</b>  | 8            | 30          | 0.038"    |

Multiloop

| One side open     | Both sides open    | Size (Fr/CH) | Length (cm)     | Guidewire |
|-------------------|--------------------|--------------|-----------------|-----------|
| <b>SS-47ML+ML</b> | <b>SOT-47ML+ML</b> | 4.7          | Multiloop 24-30 | 0.032"    |
| <b>SS-5ML+ML</b>  | <b>SOT-5ML+ML</b>  | 5            | Multiloop 24-30 | 0.035"    |
| <b>SS-6ML+ML</b>  | <b>SOT-6ML+ML</b>  | 6            | Multiloop 24-30 | 0.035"    |
| <b>SS-7ML+ML</b>  | <b>SOT-7ML+ML</b>  | 7            | Multiloop 24-30 | 0.038"    |
| <b>SS-8ML+ML</b>  | <b>SOT-8ML+ML</b>  | 8            | Multiloop 24-30 | 0.038"    |

other sizes available

Set with stainless steel guidewire add +G (for example SS-4716+G).  
Set with nitinol guidewire add +GTHS (for example SS-4716+GTHS).

## Ureter Stents

- Longlife Stent 180
- Longlife Stent 360

**Longlife Stent 180**  
Carbothane®  
Laying time up to 180 days  
Made of Polyurethane  
Tecoflex® USA

Stent getting softer  
once in the body

**Longlife Stent 360**  
100% medical Silicone  
Laying time up to 360 days  
Tecoflex® USA

Set consists of stent, pusher, suture and clamp.

### Longlife Stent 180

| Art. No.            | Size (Fr/CH) | Length (cm) | Guidewire |
|---------------------|--------------|-------------|-----------|
| <b>SOT-PRO-4826</b> | 4.8          | 26          | 0.032"    |
| <b>SOT-PRO-4828</b> | 4.8          | 28          | 0.032"    |
| <b>SOT-PRO-626</b>  | 6            | 26          | 0.035"    |
| <b>SOT-PRO-628</b>  | 6            | 28          | 0.035"    |
| <b>SOT-PRO-726</b>  | 7            | 26          | 0.038"    |
| <b>SOT-PRO-728</b>  | 7            | 28          | 0.038"    |

One side open and other sizes available

### Longlife Stent 360

|                  |   |    |        |
|------------------|---|----|--------|
| <b>SOT-S-526</b> | 5 | 26 | 0.032" |
| <b>SOT-S-528</b> | 5 | 28 | 0.032" |
| <b>SOT-S-626</b> | 6 | 26 | 0.035" |
| <b>SOT-S-628</b> | 6 | 28 | 0.035" |
| <b>SOT-S-726</b> | 7 | 26 | 0.038" |
| <b>SOT-S-728</b> | 7 | 28 | 0.038" |

One side open and other sizes available

Set with stainless steel guidewire add +G (for example SOT-S-526+G).  
Set with nitinol guidewire add +GTHS (for example SOT-S-526+GTHS).

Class IIb device – for single use only

## Ureter Stents

- Groove Stents
- Hydrophilic coated Stents

Super Drainage  
Improved flowing rate

Made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days

Groove Stent

Hydrophilic coated Stent

Groove shape



Set consists of stent, pusher, suture and clamp.

### Groove Stents

| Art. No.          | Size (Fr/CH) | Length (cm) | Guidewire |
|-------------------|--------------|-------------|-----------|
| <b>SOT-GR-726</b> | 7            | 26          | 0.038"    |
| <b>SOT-GR-728</b> | 7            | 28          | 0.038"    |

Other sizes available

### Hydrophilic coated Stents

|                   |     |    |        |
|-------------------|-----|----|--------|
| <b>SOT-H-4826</b> | 4.8 | 26 | 0.032" |
| <b>SOT-H-4828</b> | 4.8 | 28 | 0.032" |
| <b>SOT-H-626</b>  | 6   | 26 | 0.035" |
| <b>SOT-H-628</b>  | 6   | 28 | 0.035" |
| <b>SOT-H-726</b>  | 7   | 26 | 0.038" |
| <b>SOT-H-728</b>  | 7   | 28 | 0.038" |

One side open and other sizes available

Set with stainless steel guidewire add +G (for example SOT-H-626+G).  
Set with nitinol guidewire add +GTHS (for example SOT-H-626+GTHS).

Class IIb device – for single use only




Ureter Stents

- Endopyelotomy Stents
- Tumor Stents

MRI/EMG compatible

Made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days



Set consists of stent, pusher, suture and clamp.  
Stent both sides open

Endopyelotomy Stents

| Art. No.  | Size (Fr/CH) | Length (cm) | Guidewire |
|-----------|--------------|-------------|-----------|
| EP-612-24 | 6 – 12       | 24          | 0.035"    |
| EP-612-26 | 6 – 12       | 26          | 0.035"    |
| EP-612-28 | 6 – 12       | 28          | 0.035"    |
| EP-714-24 | 7 – 14       | 24          | 0.038"    |
| EP-714-26 | 7 – 14       | 26          | 0.038"    |
| EP-714-28 | 7 – 14       | 28          | 0.038"    |

Tumor Stents

|         |   |    |        |
|---------|---|----|--------|
| TSS-624 | 6 | 24 | 0.035" |
| TSS-626 | 6 | 26 | 0.035" |
| TSS-628 | 6 | 28 | 0.035" |
| TSS-724 | 7 | 24 | 0.038" |
| TSS-726 | 7 | 26 | 0.038" |
| TSS-728 | 7 | 28 | 0.038" |


Set with stainless steel guidewire add +G (for example EP-612-24+G).  
Set with nitinol guidewire add +GTHS (for example EP-612-24+GTHS).

Class IIb device – for single use only

Ureter Stents

- Tumor Stents Polyurethane
- Tumor Stents 100% Silicone

Reinforced with metal spring  
Not MRI/EMG compatible




Tumor Stent Polyurethane  
Tecoflex® USA  
Laying time up to 30 days

Reinforced with metal spring  
Not MRI/EMG compatible



Tumor Stent 100% Silicone  
Laying time up to 360 days



Set consists of stent, pusher, suture and clamp.  
Stents both sides open

Tumor Stents Polyurethane

| Art. No. | Size (Fr/CH) | Length (cm) | Guidewire |
|----------|--------------|-------------|-----------|
| TS-624   | 6            | 24          | 0.035"    |
| TS-626   | 6            | 26          | 0.035"    |
| TS-628   | 6            | 28          | 0.035"    |
| TS-724   | 7            | 24          | 0.038"    |
| TS-726   | 7            | 26          | 0.038"    |
| TS-728   | 7            | 28          | 0.038"    |

Tumor Stents 100% Silicone

|         |   |    |        |
|---------|---|----|--------|
| STS-624 | 6 | 24 | 0.035" |
| STS-626 | 6 | 26 | 0.035" |
| STS-628 | 6 | 28 | 0.035" |
| STS-724 | 7 | 24 | 0.038" |
| STS-726 | 7 | 26 | 0.038" |
| STS-728 | 7 | 28 | 0.038" |


Set with stainless steel guidewire add +G (for example TS-624+G).  
Set with nitinol guidewire add +GTHS (for example TS-624+GTHS).

Class IIb device – for single use only

Ureter Stents

- Paediatric
- Renal Transplant
- Intra Operative

Made of Polyurethane  
Tecoflex® USA  
Laying time up to  
30 days



Ureter Stents

Set consists of stent, pusher, suture and clamp. Made of Polyurethane Tecoflex® USA

Paediatric Stents

| One side open | Both sides open | Size (Fr/CH) | Length (cm) | Guidewire |
|---------------|-----------------|--------------|-------------|-----------|
| PS-310        | POS-310         | 3            | 10          | 0.018"    |
| PS-312        | POS-312         | 3            | 12          | 0.018"    |
| PS-314        | POS-314         | 3            | 14          | 0.018"    |
| PS-316        | POS-316         | 3            | 16          | 0.018"    |
| PS-410        | POS-410         | 4            | 10          | 0.025"    |
| PS-412        | POS-412         | 4            | 12          | 0.025"    |
| PS-414        | POS-414         | 4            | 14          | 0.025"    |
| PS-416        | POS-416         | 4            | 16          | 0.025"    |

Set with stainless steel guidewire add +G (for example PS-310+G).

Renal Transplant Stents

|  | Both sides open | Size (Fr/CH) | Length (cm) | Guidewire |
|--|-----------------|--------------|-------------|-----------|
|  | RS-612          | 6            | 12          | 0.035"    |
|  | RS-614          | 6            | 14          | 0.035"    |
|  | RS-616          | 6            | 16          | 0.035"    |
|  | RS-712          | 7            | 12          | 0.038"    |
|  | RS-714          | 7            | 14          | 0.038"    |
|  | RS-716          | 7            | 16          | 0.038"    |
|  | RS-812          | 8            | 12          | 0.038"    |
|  | RS-814          | 8            | 14          | 0.038"    |
|  | RS-816          | 8            | 16          | 0.038"    |

Set with stainless steel guidewire add +G (for example RS-612+G).

Intra Operative Stents both sides closed

| Size Fr/CH4.7 | Fr/CH 5 | Fr/CH 6 | Fr/CH 7 | Fr/CH 8 | Length (cm) |
|---------------|---------|---------|---------|---------|-------------|
| SSI-4716      | SSI-516 | SSI-616 | SSI-716 | SSI-816 | 16          |
| SSI-4724      | SSI-524 | SSI-624 | SSI-724 | SSI-824 | 24          |
| SSI-4726      | SSI-526 | SSI-626 | SSI-726 | SSI-826 | 26          |
| SSI-4728      | SSI-528 | SSI-628 | SSI-728 | SSI-828 | 28          |
| SSI-4730      | SSI-530 | SSI-630 | SSI-730 | SSI-830 | 30          |

Set with stainless steel guidewire

Class IIb device – for single use only

Ureter Stents

- Mono-J Stents
- Blue Stents via Ureterorenoscope

Mono-J Luer-Lock connector



Made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days

Set consists of stent, pusher, suture and clamp.

Mono-J Stents

| One side open | Both sides open | (Fr/CH) | (cm) | Guidewire |
|---------------|-----------------|---------|------|-----------|
| MJC-4770      | MJO-4770        | 4.7     | 70   | 0.032"    |
| MJC-670       | MJO-670         | 6       | 70   | 0.035"    |
| MJC-770       | MJO-770         | 7       | 70   | 0.038"    |
| MJC-870       | MJO-870         | 8       | 70   | 0.038"    |

Blue Stents via Ureterorenoscope

| One side open | Both sides open | (Fr/CH) | (cm)  | Guidewire |
|---------------|-----------------|---------|-------|-----------|
| SSU-4716      | SOU-4716        | 4.7     | 16    | 0.032"    |
| SSU-4724      | SOU-4724        | 4.7     | 24    | 0.032"    |
| SSU-4726      | SOU-4726        | 4.7     | 26    | 0.032"    |
| SSU-4728      | SOU-4728        | 4.7     | 28    | 0.032"    |
| SSU-4730      | SOU-4730        | 4.7     | 30    | 0.032"    |
| SSU-47ML+ML   | SOU-47ML+ML     | 4.7     | 24-32 | 0.032"    |

With 70cm Pusher


Set with stainless steel guidewire add +G (for example SSU-4716+G).  
Set with nitinol guidewire add +GTHS (for example SSU-4716+GTHS).

Class IIb device – for single use only

Steerable Pusher

Steerable upwards,  
downwards and turnable 360°

Easy connection



| Art. No. | Size (Fr/CH) | Steerable System for |
|----------|--------------|----------------------|
| LSDS 6   | 6            | Double J Stents      |
| LSDS 7   | 7            | Double J Stents      |


Class IIa device – for single use only



Grasping Forceps  
flexible

3-Prong, single use

To remove implanted stents or foreign bodies  
Raw material: Nitinol and Nylon66




| Art. No.   | Prongs | Size (Fr/CH) | Length (cm) |
|------------|--------|--------------|-------------|
| GF-3-3-70  | 3      | 3            | 70          |
| GF-3-3-90  | 3      | 3            | 90          |
| GF-3-3-120 | 3      | 3            | 120         |
| GF-3-4-70  | 3      | 4            | 70          |
| GF-3-4-90  | 3      | 4            | 90          |
| GF-3-4-120 | 3      | 4            | 120         |

Class IIa device

Grasping Forceps  
metallic-stiff

2-Prong, reusable

To remove implanted stents or foreign bodies  
Made of stainless steel SS304



| Art. No. | Prongs | Size (Fr/CH) | Length (cm) |
|----------|--------|--------------|-------------|
| GF 2-40  | 2      | 6            | 40          |
| GF 2-60  | 2      | 6            | 60          |

reusable-autoclavable 134°C


Class IIa device



Stainless Steel Guidewires

**PTFE coated**  
one tip stiff - one tip flexible  
Colour: PTFE dark green  
Made of stainless steel SS304

| Art. No.   | Length (cm) | Size   | Tip      |
|------------|-------------|--------|----------|
| SFT 150-32 | 150         | 0.032" | straight |
| SFT 150-35 | 150         | 0.035" | straight |
| SFT 150-38 | 150         | 0.038" | straight |
| JFT 150-32 | 150         | 0.032" | angled   |
| JFT 150-35 | 150         | 0.035" | angled   |
| JFT 150-38 | 150         | 0.038" | angled   |




other sizes available

Nitinol Super Guidewires


**75cm hydrophilic coated**  
one tip stiff - one tip flexible  
Colour: black  
Made of Nitinol

| Art. No.    | Length (cm) | Size   | Tip      |
|-------------|-------------|--------|----------|
| GTHS 150-32 | 150         | 0.032" | straight |
| GTHS 150-35 | 150         | 0.035" | straight |
| GTHA 150-32 | 150         | 0.032" | angled   |
| GTHA 150-35 | 150         | 0.035" | angled   |



Cobra Nitinol Guidewires

**75cm hydrophilic coated**  
one tip stiff - one tip flexible  
Colour: black and yellow stripes  
Made of Nitinol




Super flexible tip

| Art. No.   | Length (cm) | Size   | Tip      |
|------------|-------------|--------|----------|
| CKS 150-32 | 150         | 0.032" | straight |
| CKS 150-35 | 150         | 0.035" | straight |
| CKS 150-38 | 150         | 0.038" | straight |
| CKA 150-32 | 150         | 0.032" | angled   |
| CKA 150-35 | 150         | 0.035" | angled   |
| CKA 150-38 | 150         | 0.038" | angled   |

Ureter Catheter

To inject contrast liquid

With stylet  
Made of PEBAX  
Tecoflex® USA



Straight tip open

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| UC-O-3   | 3            | 70          |
| UC-O-4   | 4            | 70          |
| UC-O-5   | 5            | 70          |
| UC-O-6   | 6            | 70          |
| UC-O-7   | 7            | 70          |

Straight tip closed


| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| UC-R-3   | 3            | 70          |
| UC-R-4   | 4            | 70          |
| UC-R-5   | 5            | 70          |
| UC-R-6   | 6            | 70          |
| UC-R-7   | 7            | 70          |

Class IIa device – for single use only

Ureter Catheter

To inject contrast liquid

With stylet  
Made of PEBAX  
Tecoflex® USA



Angled tip open (Tiemann)

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| UC-A0-3  | 3            | 70          |
| UC-A0-4  | 4            | 70          |
| UC-A0-5  | 5            | 70          |
| UC-A0-6  | 6            | 70          |
| UC-A0-7  | 7            | 70          |

Angled tip closed (Tiemann)

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| UC-AR-3  | 3            | 70          |
| UC-AR-4  | 4            | 70          |
| UC-AR-5  | 5            | 70          |
| UC-AR-6  | 6            | 70          |
| UC-AR-7  | 7            | 70          |

Cone tip open


| Art. No. | Shaft (Fr/CH) | Cone (Fr/CH) | Length (cm) |
|----------|---------------|--------------|-------------|
| UCC-5/8  | 5             | 8            | 70          |
| UCC-6/10 | 6             | 10           | 70          |
| UCC-7/12 | 7             | 12           | 70          |

Class IIa device – for single use only

# Ureter Access Sheath Standard

Made of Polyurethane  
Tecoflex® USA  
Dilator: PEBAX

Provides **ureteral dilation** and  
a **continous working channel**  
for the introduction of Endoscopes  
and Instruments during ureteral  
access procedures.




| Art. No.  | Sheath Inner Diameter<br>Size (Fr/CH) = mm |   |     | Sheath Length |
|-----------|--------------------------------------------|---|-----|---------------|
| UAS 10-13 | 10                                         | = | 3.3 | 13            |
| UAS 10-20 | 10                                         | = | 3.3 | 20            |
| UAS 10-28 | 10                                         | = | 3.3 | 28            |
| UAS 10-30 | 10                                         | = | 3.3 | 30            |
| UAS 10-35 | 10                                         | = | 3.3 | 35            |
| UAS 10-45 | 10                                         | = | 3.3 | 45            |
| UAS 12-13 | 12                                         | = | 4   | 13            |
| UAS 12-20 | 12                                         | = | 4   | 20            |
| UAS 12-28 | 12                                         | = | 4   | 28            |
| UAS 12-30 | 12                                         | = | 4   | 30            |
| UAS 12-35 | 12                                         | = | 4   | 35            |
| UAS 12-45 | 12                                         | = | 4   | 45            |
| UAS 14-13 | 14                                         | = | 4.7 | 13            |
| UAS 14-20 | 14                                         | = | 4.7 | 20            |
| UAS 14-28 | 14                                         | = | 4.7 | 28            |
| UAS 14-30 | 14                                         | = | 4.7 | 30            |
| UAS 14-35 | 14                                         | = | 4.7 | 35            |
| UAS 14-45 | 14                                         | = | 4.7 | 45            |

Class IIa device – for single use only

# Ureter Access Sheath Braided

Made of Polyurethane (braided)  
Tecoflex® USA  
Dilator: PEBAX  
Silver part: Stainless steel, braided tube

Provides **ureteral dilation** and  
a **continous working channel**  
for the introduction of Endoscopes  
and Instruments during ureteral  
access procedures.



| Art. No.    | Sheath Inner Diameter<br>Size (Fr/CH) = mm |   |     | Sheath Length |
|-------------|--------------------------------------------|---|-----|---------------|
| UAS-B 10-13 | 10                                         | = | 3.3 | 13            |
| UAS-B 10-20 | 10                                         | = | 3.3 | 20            |
| UAS-B 10-28 | 10                                         | = | 3.3 | 28            |
| UAS-B 10-30 | 10                                         | = | 3.3 | 30            |
| UAS-B 10-35 | 10                                         | = | 3.3 | 35            |
| UAS-B 10-45 | 10                                         | = | 3.3 | 45            |
| UAS-B 12-13 | 12                                         | = | 4   | 13            |
| UAS-B 12-20 | 12                                         | = | 4   | 20            |
| UAS-B 12-28 | 12                                         | = | 4   | 28            |
| UAS-B 12-30 | 12                                         | = | 4   | 30            |
| UAS-B 12-35 | 12                                         | = | 4   | 35            |
| UAS-B 12-45 | 12                                         | = | 4   | 45            |
| UAS-B 14-13 | 14                                         | = | 4.7 | 13            |
| UAS-B 14-20 | 14                                         | = | 4.7 | 20            |
| UAS-B 14-28 | 14                                         | = | 4.7 | 28            |
| UAS-B 14-30 | 14                                         | = | 4.7 | 30            |
| UAS-B 14-35 | 14                                         | = | 4.7 | 35            |
| UAS-B 14-45 | 14                                         | = | 4.7 | 45            |


UAS-BH hydrophilic coated

Class IIa device – for single use only




Kidney Stone Basket  
Standard


Tube made of PTFE  
Basket and wire made of Nitinol



**Helical**  
Add to Art. No. -H  
(Example SE 3-470-H)



**Flatwire**  
Add to Art. No. -F  
(Example SE 3-470-F)



**Tipless**  
Add to Art. No. -T  
(Example SE 3-470-T)

| Art. No.  | Size (Fr/CH) | Number of wires | Length (cm) |
|-----------|--------------|-----------------|-------------|
| SE 3-470  | 3            | 4               | 70          |
| SE 3-490  | 3            | 4               | 90          |
| SE 3-4120 | 3            | 4               | 120         |
| SE 4-470  | 4            | 4               | 70          |
| SE 4-490  | 4            | 4               | 90          |
| SE 4-4120 | 4            | 4               | 120         |

For flexible Ureterorenoscope

|                    |     |   |     |
|--------------------|-----|---|-----|
| NITI-DIS 1.8-4-120 | 1.8 | 4 | 120 |
| NITI-DIS 2.4-4-120 | 2.4 | 4 | 120 |

Baskets tipless

Class IIa device – for single use only

G-PAW Stone Basket

Nitinol

G-Paw is designed to offer advantages of both a stone grasper and a stone basket. It allows easy grasp, release and removal of stones of various sizes throughout the urinary tract.

**D-Shape Handle**  
Add to Art. No. -D  
(Example BUSB9117370-D)



**Fish Type Handle**  
Add to Art. No. -F  
(Example BUSB9117370-F)



**Scissor Type Handle**  
Add to Art. No. -S  
(Example BUSB9117370-S)



Stone Basket (Nitinol-G-Paw)

| Art. No.    | Size (Fr/CH) | Number of wires | Length (cm) |
|-------------|--------------|-----------------|-------------|
| BUSB9117370 | 1.7          | 3               | 70          |
| BUSB9119370 | 1.9          | 3               | 70          |
| BUSB9122370 | 2.2          | 3               | 70          |
| BUSB9125370 | 2.5          | 3               | 70          |
| BUSB9130370 | 3.0          | 3               | 70          |
| BUSB9117390 | 1.7          | 3               | 90          |
| BUSB9119390 | 1.9          | 3               | 90          |
| BUSB9122390 | 2.2          | 3               | 90          |
| BUSB9125390 | 2.5          | 3               | 90          |
| BUSB9130390 | 3.0          | 3               | 90          |
| BUSB9117310 | 1.7          | 3               | 120         |
| BUSB9119310 | 1.9          | 3               | 120         |
| BUSB9122310 | 2.2          | 3               | 120         |
| BUSB9125310 | 2.5          | 3               | 120         |
| BUSB9130310 | 3.0          | 3               | 120         |



# X-Circle Stone Basket

## Nitinol

X-Circle provides a unique net type geometry to capture and extract intact stones and prevent migration.

**D-Shape Handle**  
Add to Art. No. -D  
(Example BUSX2117370-D)



**Fish Type Handle**  
Add to Art. No. -F  
(Example BUSX2117370-F)



**Scissor Type Handle**  
Add to Art. No. -S  
(Example BUSX2117370-S)



### Stone Basket (Nitinol-X-Circle)

| Art. No.    | Size (Fr/CH) | Number of wires | Length (cm) |
|-------------|--------------|-----------------|-------------|
| BUSX2117370 | 1.7          | 12              | 70          |
| BUSX2119370 | 1.9          | 12              | 70          |
| BUSX2122370 | 2.2          | 12              | 70          |
| BUSX2125370 | 2.5          | 12              | 70          |
| BUSX2130370 | 3.0          | 12              | 70          |
| BUSX2117390 | 1.7          | 12              | 90          |
| BUSX2119390 | 1.9          | 12              | 90          |
| BUSX2122390 | 2.2          | 12              | 90          |
| BUSX2125390 | 2.5          | 12              | 90          |
| BUSX2130390 | 3.0          | 12              | 90          |
| BUSX2117310 | 1.7          | 12              | 120         |
| BUSX2119310 | 1.9          | 12              | 120         |
| BUSX2122310 | 2.2          | 12              | 120         |
| BUSX2125310 | 2.5          | 12              | 120         |
| BUSX2130310 | 3.0          | 12              | 120         |

# Kidney Stone Stop

## Nitinol

Stone Stop is used to prevent stone migration to the kidney during lithotripsy.

**Stone Stop**




### Stone Basket (Stone Stop)

| Art. No.    | Size (Fr/CH) | Basket Ø (mm) | Length (cm) |
|-------------|--------------|---------------|-------------|
| BUSS0107090 | 3            | 7             | 90          |
| BUSS0110090 | 3            | 10            | 90          |
| BUSS0107110 | 3            | 7             | 110         |
| BUSS0110110 | 3            | 10            | 110         |
| BUSS0107120 | 3            | 7             | 120         |
| BUSS0110120 | 3            | 10            | 120         |

# Nephrostomy Drainage Set with Pigtail Catheter

Made of Polyurethane  
Tecoflex® USA



# Nephrostomy Drainage Set with Pigtail Catheter

A convenient complete set for percutaneous placement of drainage catheter.  
Soft polyurethane material provides superior biocompatibility and drainage.  
Highly radiopaque for precise placement.  
Laying time up to 30 days

- 1 Set consists of:
- PIGTAIL – Drainage Catheter 35cm, usable length 30cm, plus straightener
  - 2-part puncture needle 1,3mm (18 GA)
  - Lunderquist Guidewire 80cm J-Tip, fixed core, stainless steel, 0.035"
  - Fascial Dilators 20cm, ABS
  - Urine Bag connector

## Set with PIGTAIL Catheter

Catheter hydrophilic coated

| Art. No. | Catheter Size (Fr/CH) | Fascial Dilators (Fr/CH) |
|----------|-----------------------|--------------------------|
| PDS-P6   | 6                     | 5 – 6 – 7                |
| PDS-P7   | 7                     | 6 – 7 – 8                |
| PDS-P8   | 8                     | 7 – 8 – 9                |
| PDS-P9   | 9                     | 7 – 8 – 9 – 10           |
| PDS-P10  | 10                    | 6 – 8 – 10 – 12          |
| PDS-P12  | 12                    | 6 – 8 – 10 – 12          |
| PDS-P14  | 14                    | 6 – 8 – 12 – 14          |

## PIGTAIL Catheter only

| Art. No. | Size (Fr/CH) |
|----------|--------------|
| PCN-P6   | 6            |
| PCN-P7   | 7            |
| PCN-P8   | 8            |
| PCN-P9   | 9            |
| PCN-P10  | 10           |
| PCN-P12  | 12           |
| PCN-P14  | 14           |


## PIGTAIL Exchange Set

| Art. No.  | Catheter Size (Fr/CH) |
|-----------|-----------------------|
| PDS-PE 8  | 8                     |
| PDS-PE 10 | 10                    |
| PDS-PE 12 | 12                    |
| PDS-PE 14 | 14                    |

Same as drainage set without puncture  
needle

# Nephrostomy Drainage Set with Malecot Catheter

Made of Polyurethane  
Tecoflex® USA



# Nephrostomy Drainage Set with Malecot Catheter

A convenient complete set for percutaneous placement of drainage catheter.  
Soft polyurethane material provides superior biocompatibility and drainage.  
Highly radiopaque for precise placement.  
Laying time up to 30 days

- 1 Set consists of:
- MALECOT – Drainage Catheter 35cm, usable length 30cm, plus straightener
  - 2-part puncture needle 1,3mm (18 GA)
  - Lunderquist Guidewire 80cm J-Tip, fixed core, stainless steel, 0.035"
  - Fascial Dilators 20cm, ABS
  - Urine Bag connector

## Set with MALECOT Catheter

Catheter hydrophilic coated

| Art. No.                      | Catheter Size (Fr/CH) | Fascial Dilators (Fr/CH) |
|-------------------------------|-----------------------|--------------------------|
| PDS-M10                       | 10                    | 6 – 8 – 10 – 12          |
| PDS-M12                       | 12                    | 6 – 8 – 12 – 14          |
| PDS-M14                       | 14                    | 8 – 12 – 14 – 16         |
| PDS-M16                       | 16                    | 8 – 12 – 14 – 16         |
| PDS-M18 / M20 / M22 available |                       |                          |


## MALECOT Catheter only

| Art. No.                      | Size (Fr/CH) |
|-------------------------------|--------------|
| PCN-M10                       | 10           |
| PCN-M12                       | 12           |
| PCN-M14                       | 14           |
| PCN-M16                       | 16           |
| PCN-M18 / M20 / M22 available |              |

Class IIb device – for single use only

# Nephrostomy Drainage Set

with 100% Silicone Integral Catheter



# Nephrostomy Drainage Set

with 100% Silicone Integral Catheter

A convenient complete set for percutaneous placement of a 2-way Silicone Integral Balloon Catheter. Excellent biocompatibility and drainage.  
Laying time up to 30 days

- 1 Set consists of:
- Drainage Catheter 40cm, usable length 35cm, plus straightener
  - 2-part puncture needle 1,3mm (18 GA)
  - Lunderquist Guidewire 80cm J-Tip, fixed core, stainless steel, 0.035"
  - Fascial Dilators 20cm, ABS
  - Urine Bag connector

## Set with 2-way INTEGRAL Balloon Catheter

| Art. No.   | Catheter Size (Fr/CH) | Balloon (cc/ml) | Facial Dilators (Fr/CH) |
|------------|-----------------------|-----------------|-------------------------|
| CH-PCNL 10 | 10                    | 1.5             | 6 – 8 – 10              |
| CH-PCNL 12 | 12                    | 2               | 8 – 10 – 12             |
| CH-PCNL 14 | 14                    | 3               | 10 – 12 – 14            |
| CH-PCNL 16 | 16                    | 5               | 10 – 12 – 14 – 16       |

## INTEGRAL Exchange Set

Same as drainage set without puncture needle

|              |    |     |                   |
|--------------|----|-----|-------------------|
| CH-PCNL 10-E | 10 | 1.5 | 6 – 8 – 10        |
| CH-PCNL 12-E | 12 | 2   | 8 – 10 – 12       |
| CH-PCNL 14-E | 14 | 3   | 10 – 12 – 14      |
| CH-PCNL 16-E | 16 | 5   | 10 – 12 – 14 – 16 |

## INTEGRAL Balloon Catheter only

|              |    |     |
|--------------|----|-----|
| CH-PCNL 10-S | 10 | 1.5 |
| CH-PCNL 12-S | 12 | 2   |
| CH-PCNL 14-S | 14 | 3   |
| CH-PCNL 16-S | 16 | 5   |
| CH-PCNL 18-S | 18 | 5   |
| CH-PCNL 20-S | 20 | 5   |


Class IIb device – for single use only



# Nephrostomy Drainage

## Pigtail Catheter with Locking Mechanism

Catheter made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days



**Set with Pigtail Catheter  
with Locking Mechanism**

| Art. No.   | Size (Fr/CH) | Length (cm) |
|------------|--------------|-------------|
| PDS-P8-LM  | 8            | 40          |
| PDS-P9-LM  | 9            | 40          |
| PDS-P10-LM | 10           | 40          |
| PDS-P12-LM | 12           | 40          |
| PDS-P14-LM | 14           | 40          |
| PDS-P16-LM | 16           | 40          |

**Pigtail Catheter  
with Locking Mechanism single**

| Art. No.   | Size (Fr/CH) | Length (cm) |
|------------|--------------|-------------|
| PCN-P6-LM  | 6            | 40          |
| PCN-P7-LM  | 7            | 40          |
| PCN-P8-LM  | 8            | 40          |
| PCN-P9-LM  | 9            | 40          |
| PCN-P10-LM | 10           | 40          |
| PCN-P12-LM | 12           | 40          |
| PCN-P14-LM | 14           | 40          |

Class IIb device – for single use only

# Nephrostomy Drainage Set

## Direct Puncture


Catheter made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days



Scalpel



Initial Puncture Needle to be assembled with Pigtail-Catheter



Pigtail – Drainage Catheter, Hydrophilic coated

**Direct Puncture Set  
with Pigtail Catheter**

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| DPD-7    | 7            | 35          |
| DPD-8    | 8            | 35          |
| DPD-9    | 9            | 35          |
| DPD-10   | 10           | 35          |

**Direct Puncture Set  
with Pigtail Catheter  
with Locking Mechanism**

| Art. No.  | Size (Fr/CH) | Length (cm) |
|-----------|--------------|-------------|
| DPD-8-LM  | 8            | 40          |
| DPD-9-LM  | 9            | 40          |
| DPD-10-LM | 10           | 40          |
| DPD-12-LM | 12           | 40          |
| DPD-14-LM | 14           | 40          |
| DPD-16-LM | 16           | 40          |

Class IIb device – for single use only

# Re-Entry MALECOT Drainage Catheter

Used after PCNL procedure  
For Renal Drainage  
With Urine Bag connector

Made of Polyurethane  
Tecoflex® USA  
Laying time up to 30 days



| Art. No. | Size (Fr/CH) |
|----------|--------------|
| RM 10/6  | 10           |
| RM 12/6  | 12           |
| RM 14/6  | 14           |
| RM 16/8  | 16           |
| RM 18/8  | 18           |

Class IIb device – for single use only

# Initial Puncture Needles


Used for percutaneous Nephrostomy and all other percutaneous applications  
where a 0.035" guidewire has to be used.

## 2-Part Puncture Needle, Trocar tip



| Art. No.              | Size (GA) | Length (cm) |
|-----------------------|-----------|-------------|
| IPN-18-20-2T          | 18        | 20          |
| other sizes available |           |             |

## 3-Part Puncture Needle, Bevel tip



| Art. No.              | Size (GA) | Length (cm) |
|-----------------------|-----------|-------------|
| IPN-18-20-3B          | 18        | 20          |
| other sizes available |           |             |

# Chiba Needles



| Art. No.              | Size (GA) | Length (cm) |
|-----------------------|-----------|-------------|
| CH-22-20              | 22        | 20          |
| CH-22-28              | 22        | 28          |
| other sizes available |           |             |

Class IIa device – for single use only

Fascial Dilators

Hydrophilic coated

With female luer-lock  
Accepts 0.038" guidewire

Set

| Art. No. | Dilators | Size (Fr/CH)       | Length (cm) |
|----------|----------|--------------------|-------------|
| FDS 6-10 | 5        | 6 – 7 – 8 – 9 – 10 | 20          |


|                       |   |                       |    |
|-----------------------|---|-----------------------|----|
| FDS 8-16              | 5 | 8 – 10 – 12 – 14 – 16 | 20 |
| other sizes available |   |                       |    |

Single 20cm

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| FD 4-20  | 4            | 20          |
| FD 6-20  | 6            | 20          |
| FD 8-20  | 8            | 20          |
| FD 10-20 | 10           | 20          |
| FD 12-20 | 12           | 20          |
| FD 14-20 | 14           | 20          |
| FD 16-20 | 16           | 20          |
| FD 18-20 | 18           | 20          |
| FD 20-20 | 20           | 20          |

Single 30cm


|       |    |    |
|-------|----|----|
| FD 8  | 8  | 30 |
| FD 10 | 10 | 30 |
| FD 12 | 12 | 30 |
| FD 14 | 14 | 30 |
| FD 16 | 16 | 30 |



Ureter Dilator Set

70cm length

Hydrophilic coated



Art. No.  
UDS 6/16-70

Ureter Balloon-Dilator



| Art. No. | Shaft size (Fr/CH) | Total length | Balloon diam. | Balloon length |
|----------|--------------------|--------------|---------------|----------------|
| BD 5/4   | 5                  | 70cm         | CH 15/5mm     | 40mm           |

Accepts 0.035" guidewire, Supplied with stopcock, 18atm rated burst pressure

Transureteroscopic Balloon-Dilator

| Art. No. | Shaft size (Fr/CH) | Total length | Balloon diam. | Balloon length |
|----------|--------------------|--------------|---------------|----------------|
| TU 4/4   | 4                  | 70cm         | CH 12/3.5mm   | 40mm           |

Provided with stylet, Accepts 0.018" guidewire, Supplied with stopcock  
18atm rated burst pressure

Nephrostomy Track Balloon-Dilator



| Art. No.    | Balloon diam. | Balloon length | Total length |                     |
|-------------|---------------|----------------|--------------|---------------------|
| NT 10/12-AM | 30Fr/CH       | 120mm          | 550mm        | with Amplatz Sheath |

Accepts 0.035" guidewire, 12atm rated burst pressure  
Supplied with Stopcock and Amplatz Sheath

Inflation Device

for Inflation of urological balloon dilators  
20cc capacity  
Includes 15cm Extension tube  
With stop-cock and luer-lock




|          |
|----------|
| Art. No. |
| ID 20    |



# Nottingham One-Step Dilator

Hydrophilic coated, to dilate ureter


Total length 70cm, Taper length 6cm, accepts 0.038" guidewire



| Art. No. | Shaft Size (Fr/CH) | Taper Size (Fr/CH) from |
|----------|--------------------|-------------------------|
| NOT 6-10 | 10                 | 6 to 10                 |
| NOT 6-12 | 12                 | 6 to 12                 |

# Female Dilator

To dilate the female urethra




| Art. No. | Taper Size (Fr/CH) from | Length |
|----------|-------------------------|--------|
| FUD      | 3 to 10mm               | 100mm  |

# Meatal Dilator

To dilate urethral meatus

Paediatric or adult model available



| Art. No. | Taper Size (Fr/CH) from | Length |
|----------|-------------------------|--------|
| MDA      | 2 to 8mm                | 45mm   |
| MDP      | 1 to 4mm                | 20mm   |

# Urethra Dilator Set curved-soft

For the dilation of the urethra

Set with 3 dilators  
12Fr, 16Fr and 20Fr

Available with  
open or closed tip




Soft dilators  
Made of Polyurethane  
Tecoflex® USA

| Art. No. | Tip    | Length (cm) |
|----------|--------|-------------|
| URS-C-40 | closed | 40          |
| URS-O-40 | open   | 40          |

# Amplatz Renal Dilation Set

## Hydrophilic coated

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.



| Art. No.  | Size (Fr/CH) |
|-----------|--------------|
| ARDS 8/30 | 8 to 30      |

Guiding Catheter 8Fr/CH  
Length 80cm

# Amplatz Renal Dilators

## Hydrophilic coated

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.

### Single

| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| RD 18    | 18           | 30          |
| RD 20    | 20           | 30          |
| RD 22    | 22           | 30          |
| RD 24    | 24           | 30          |
| RD 26    | 26           | 30          |
| RD 28    | 28           | 30          |
| RD 30    | 30           | 30          |

# Amplatz Sheath

## Hydrophilic coated

Used during renal dilatation to provide atraumatic working tracks after removal of the dilator.



| Art. No. | Size (Fr/CH) | Length (cm) |
|----------|--------------|-------------|
| AM 10    | 10           | 30          |
| AM 24    | 24           | 17          |
| AM 26    | 26           | 17          |
| AM 28    | 28           | 17          |
| AM 30    | 30           | 17          |

# Mini PCNL – Dilation Set

## Percutaneous Nephrolithotomy

Dilators made of Teflon/PTFE

### Amplatz Sheath

12Fr, Ø ID 4.1mm  
14Fr, Ø ID 4.7mm  
16Fr, Ø ID 5.4mm

Set consists of  
1 Screw Dilator  
12Fr, 14Fr or 16Fr  
Length 30cm  
1 Stainless Steel Guidewire  
“super-stiff” 80cm, 0,038”  
1 2-Part Needle  
1 Amplatz Sheath  
12Fr, 14Fr or 16Fr  
Length 17cm

guidewire PTFE  
coated

Art. No.  
**PCNL-M 12**  
**PCNL-M 14**  
**PCNL-M 16**

Class IIa device – for single use only

# Semi-Automatic Biopsy-Gun

## Gun and Needle combined

To take a 10mm or 20mm  
biopsy sample from Prostate

To be seen under fluoroscopy  
or ultrasound



| Art. No.         | Size (GA) | Length (cm) |
|------------------|-----------|-------------|
| <b>SEMI 1410</b> | 14        | 10          |
| <b>SEMI 1415</b> | 14        | 15          |
| <b>SEMI 1420</b> | 14        | 20          |
| <b>SEMI 1610</b> | 16        | 10          |
| <b>SEMI 1615</b> | 16        | 15          |
| <b>SEMI 1620</b> | 16        | 20          |
| <b>SEMI 1810</b> | 18        | 10          |
| <b>SEMI 1815</b> | 18        | 15          |
| <b>SEMI 1820</b> | 18        | 20          |


other sizes available



Pos. 0: Initial position      Pos. 1: for 10mm sample      Pos. 2: for 20mm sample

Class IIa device – for single use only

Suprapubic Drainage  
Cystocath



Sharp metal tip

Suprapubic Drainage  
Cystocath

Set with Balloon Catheter 40cm

Set consists of: 2-way Suprapubic Catheter, Scalpel, Trocar and Sleeve

| Art. No. | Size (Fr/CH) | Catheter Length (cm) |
|----------|--------------|----------------------|
| SDS 10   | 10           | 40                   |
| SDS 12   | 12           | 40                   |
| SDS 14   | 14           | 40                   |
| SDS 16   | 16           | 40                   |
| SDS 18   | 18           | 40                   |

Puncture Set without Balloon Catheter


Set consists of: Scalpel, Trocar and Sleeve

| Art. No. | Size (Fr/CH) |
|----------|--------------|
| SPS 10   | 10           |
| SPS 12   | 12           |
| SPS 14   | 14           |
| SPS 16   | 16           |
| SPS 18   | 18           |



# Suprapubic Drainage Integral

Accepts 0.038" guidewire  
100% Silicone, transparent



**Puncture Set INTEGRAL**

| Art. No.         | Size (Fr/CH) | Catheter Length (cm) |
|------------------|--------------|----------------------|
| <b>SDS-MB 12</b> | 12           | 40                   |
| <b>SDS-MB 14</b> | 14           | 40                   |

Set consists of 2-way Suprapubic Catheter with integrated Balloon and Metal Trocar

**Exchange Set INTEGRAL**


| Art. No.         | Size (Fr/CH) | Catheter Length (cm) |
|------------------|--------------|----------------------|
| <b>SDS-BX 14</b> | 14           | 40                   |

Set consists of 2-way Suprapubic Catheter with integrated Balloon and guidewire

# Suprapubic Drainage with Pigtail Catheter

Set consists of:

- Trocar
- PTFE coated stainless steel guidewire 0.038", 80cm
- Catheter-Straightener stainless steel
- Pigtail Catheter Polyurethane Tecoflex®



**Drainage Set with Pigtail Catheter**

| Art. No.         | Size (Fr/CH) |
|------------------|--------------|
| <b>SDS-MP 10</b> | 10           |
| <b>SDS-MP 12</b> | 12           |
| <b>SDS-MP 14</b> | 14           |

# Haematuria Catheter


Silicone Elastomer, blue

2- and 3-way for Prostate operation

Shaft reinforced with Nylon Suture to prevent collapsing

Catheter 40cm

Balloon 60ml/cc



# Haematuria Catheter

Silicone Elastomer, blue

2-way NELATON tip

| Art. No.  | Size (Fr/CH) | Balloon (ml) |
|-----------|--------------|--------------|
| 9387-1860 | 18           | 60           |
| 9387-2060 | 20           | 60           |
| 9387-2260 | 22           | 60           |
| 9387-2460 | 24           | 60           |

3-way NELATON tip

| Art. No.  | Size (Fr/CH) | Balloon (ml) |
|-----------|--------------|--------------|
| 9390-1860 | 18           | 60           |
| 9390-2060 | 20           | 60           |
| 9390-2260 | 22           | 60           |
| 9390-2460 | 24           | 60           |

2-way COUVELAIRE tip

|           |    |    |
|-----------|----|----|
| 9394-1860 | 18 | 60 |
| 9394-2060 | 20 | 60 |
| 9394-2260 | 22 | 60 |
| 9394-2460 | 24 | 60 |

3-way COUVELAIRE tip

|           |    |    |
|-----------|----|----|
| 9395-1860 | 18 | 60 |
| 9395-2060 | 20 | 60 |
| 9395-2260 | 22 | 60 |
| 9395-2460 | 24 | 60 |

2-way DUFOUR tip

|           |    |    |
|-----------|----|----|
| 9391-1860 | 18 | 60 |
| 9391-2060 | 20 | 60 |
| 9391-2260 | 22 | 60 |
| 9391-2460 | 24 | 60 |

3-way DUFOUR tip

|           |    |    |
|-----------|----|----|
| 9396-1860 | 18 | 60 |
| 9396-2060 | 20 | 60 |
| 9396-2260 | 22 | 60 |
| 9396-2460 | 24 | 60 |

# Haematuria Catheter

100% Silicone, transparent

3-way Nelaton tip

| Art. No.    | Size (Fr/CH) | Balloon (ml) |
|-------------|--------------|--------------|
| 9390VS-1630 | 16           | 30           |
| 9390VS-1830 | 18           | 30           |
| 9390VS-2030 | 20           | 30           |
| 9390VS-2230 | 22           | 30           |
| 9390VS-2430 | 24           | 30           |

TUR Electrodes

For Prostate Surgery



Sphere Electrode (5mm Ball)



Cutting Loop for 30° Optic



Necklace Electrode



Cold Knife Electrode



Roller Electrode



Kriloktrode Electrode  
(Simultaneous Cutting & Vaporization)



Collins Knife



Spike Roller Electrode

TUR Electrodes MONOPOLAR

STORZ, 2 Stem, 30° Degrees, 24Fr

|                    |           |
|--------------------|-----------|
| Cutting Loop       | W-DC-24Fr |
| Collin's Knife     | W-DZ-24Fr |
| Roller Electrode   | W-DR-24Fr |
| Ball Electrode 3mm | W-DB-24Fr |
| Ball Electrode 5mm | W-DS-24Fr |
| Kriloktrode        | W-DK-24Fr |
| Cold Knife         | W-DP-24Fr |
| Necklace           | W-DN-24Fr |

|           |           |
|-----------|-----------|
| Spike 3mm | W-DW-24Fr |
| Spike 5mm | W-DV-24Fr |

OLYMPUS, 1 Stem, 12° Degrees, 24Fr

|                  |           |
|------------------|-----------|
| Cutting Loop     | W-OC-24Fr |
| Collin's Knife   | W-OZ-24Fr |
| Roller Electrode | W-OR-24Fr |

|                      |           |
|----------------------|-----------|
| Ball Electrode 3mm   | W-OB-24Fr |
| Sphere Electrode 5mm | W-OS-24Fr |
| Kriloktrode          | W-OK-24Fr |

WOLF, 1 Stem, 15° Degrees, 24Fr

|                  |           |
|------------------|-----------|
| Cutting Loop     | W-WC-24Fr |
| Collin's Knife   | W-WZ-24Fr |
| Roller Electrode | W-WR-24Fr |

|                      |           |
|----------------------|-----------|
| Ball Electrode 3mm   | W-WB-24Fr |
| Sphere Electrode 5mm | W-WS-24Fr |
| Kriloktrode          | W-WK-24Fr |

STORZ, 1 Stem, 30° Degrees, 24Fr

|                    |           |
|--------------------|-----------|
| Cutting Loop       | W-SC-24Fr |
| Collin's Knife     | W-SZ-24Fr |
| Roller Electrode   | W-SR-24Fr |
| Ball Electrode 3mm | W-SB-24Fr |
| Ball Electrode 5mm | W-SS-24Fr |
| Kriloktrode        | W-SK-24Fr |
| Cold Knife         | W-SP-24Fr |
| Necklace           | W-SN-24Fr |

|           |           |
|-----------|-----------|
| Spike 3mm | W-SW-24Fr |
| Spike 5mm | W-SV-24Fr |

OLYMPUS, 1 Stem, 30° Degrees, 24Fr

|                  |           |
|------------------|-----------|
| Cutting Loop     | W-LC-24Fr |
| Collin's Knife   | W-LZ-24Fr |
| Roller Electrode | W-LR-24Fr |

|                      |           |
|----------------------|-----------|
| Ball Electrode 3mm   | W-LB-24Fr |
| Sphere Electrode 5mm | W-LS-24Fr |
| Kriloktrode          | W-LK-24Fr |

ACMI

|              |        |
|--------------|--------|
| Cutting Loop | UAC 24 |
| Cutting Loop | UAC 26 |
| Cutting Loop | UAC 28 |

TUR Electrodes BIPOLAR

|                                     |                       |                                     |                   |
|-------------------------------------|-----------------------|-------------------------------------|-------------------|
| STORZ, 2 Stem, 30° Degrees, 24/26Fr |                       | STORZ, 2 Stem, 30° Degrees, 24/26Fr |                   |
| Cutting Loop                        | UBSCP                 | Coagulation Electrode               | UBSCL (pointed)   |
| Cutting Loop                        | UBSCP40 (0.40mm wire) | Coagulation Electrode               | UBSCLV (pointed)  |
| Cutting Loop                        | UBSC                  | Half moon Vaporization              |                   |
| Cutting Loop                        | UBSC30 (0.30mm wire)  | Electrode, ball-shaped              | UBSB              |
| Cutting Loop                        | UBSC40 (0.40mm wire)  |                                     |                   |
| Cutting Loop                        | UBSCDV (small)        |                                     |                   |
| Cutting Loop                        | UBSCD (small)         |                                     |                   |
| Cutting Loop                        | UBSCB (straight)      |                                     |                   |
| OLYMPUS, 12° Degrees, 24Fr          |                       | OLYMPUS, 30° Degrees, 24Fr          |                   |
| Cutting Loop                        | UBOC12-S (small)      | Cutting Loop                        | UBOC30-S (small)  |
| Cutting Loop                        | UBOC12-M (medium)     | Cutting Loop                        | UBOC30-M (medium) |
| OLYMPUS, 12° and 30° Degrees, 24Fr  |                       |                                     |                   |
| Angeled Loop                        | UBOCA-S (small)       |                                     |                   |
| Angeled Loop                        | UBOCA-M (medium)      |                                     |                   |
| OLYMPUS, 12° and 30° Degrees        |                       | OLYMPUS, 12°, 30° and 45° Degrees   |                   |
| Roller Electrode                    | UBOR                  | Needle,                             |                   |
| Loop Electrode                      | UBOC12-D, 12° large   | Angeled Loop                        |                   |
| Loop Electrode                      | UBOC30-D, 30° large   | UBOL                                |                   |
| Band Electrode                      | UBOC12-T, 12° large   |                                     |                   |
| Band Electrode                      | UBOC30-T, 30° large   |                                     |                   |
| Button for Plasma                   | UBOB                  |                                     |                   |

Bladder Evaquator

Used to perform irrigation and tissue collection during TUR prostate or bladder surgery

ELICK



Art. No.  
**EEC**  
For single use only




## Bladder Irrigation Set



Boxes of 25 pcs.

Class IIa device – for single use only

## TUR-Collector



## Urine Bag

2000 cc capacity  
with 110cm tube  
drop-chamber  
anti-bacteria plug  
cross-valve



2-Way Balloon Catheters



2-Way Balloon Catheter

Longlife Silicone Elastomer

Nelaton tip

| Art. No.   | Size<br>(Fr/CH) | Balloon<br>(cc/ml) | Length<br>(cm) |
|------------|-----------------|--------------------|----------------|
| 9385S-0603 | 6               | 3                  | 20             |
| 9385S-0805 | 8               | 5                  | 20             |
| 9385S-1005 | 10              | 5                  | 20             |
| 9385S-1210 | 12              | 5 to 15            | 40             |
| 9385S-1410 | 14              | 5 to 15            | 40             |
| 9385S-1610 | 16              | 5 to 15            | 40             |
| 9385S-1810 | 18              | 5 to 15            | 40             |
| 9385S-2010 | 20              | 5 to 15            | 40             |
| 9385S-2210 | 22              | 5 to 15            | 40             |
| 9385S-2410 | 24              | 5 to 15            | 40             |
| 9385S-2610 | 26              | 5 to 15            | 40             |

Tiemann tip

| Art. No.   | Size<br>(Fr/CH) | Balloon<br>(cc/ml) | Length<br>(cm) |
|------------|-----------------|--------------------|----------------|
| 9310S-1210 | 12              | 5 to 15            | 40             |
| 9310S-1410 | 14              | 5 to 15            | 40             |
| 9310S-1610 | 16              | 5 to 15            | 40             |
| 9310S-1810 | 18              | 5 to 15            | 40             |
| 9310S-2010 | 20              | 5 to 15            | 40             |

Suprapubic tip

|          |    |         |    |
|----------|----|---------|----|
| 9399S-12 | 12 | 5 to 15 | 40 |
| 9399S-14 | 14 | 5 to 15 | 40 |
| 9399S-16 | 16 | 5 to 15 | 40 |
| 9399S-18 | 18 | 5 to 15 | 40 |
| 9399S-20 | 20 | 5 to 15 | 40 |
| 9399S-22 | 22 | 5 to 15 | 40 |
| 9399S-24 | 24 | 5 to 15 | 40 |

100% Silicone Nelaton

Nelaton tip

|             |    |    |    |
|-------------|----|----|----|
| 9385VS-1210 | 12 | 10 | 40 |
| 9385VS-1410 | 14 | 10 | 40 |
| 9385VS-1610 | 16 | 10 | 40 |
| 9385VS-1810 | 18 | 10 | 40 |
| 9385VS-2010 | 20 | 10 | 40 |
| 9385VS-2210 | 22 | 10 | 40 |
| 9385VS-2410 | 24 | 10 | 40 |

100% Silicone Suprapubic

Suprapubic tip

|           |    |    |    |
|-----------|----|----|----|
| 9399VS-12 | 12 | 10 | 40 |
| 9399VS-14 | 14 | 10 | 40 |
| 9399VS-16 | 16 | 10 | 40 |
| 9399VS-18 | 18 | 10 | 40 |
| 9399VS-20 | 20 | 10 | 40 |
| 9399VS-22 | 22 | 10 | 40 |
| 9399VS-24 | 24 | 10 | 40 |

## Catheter Gel

**For easy insertion of a catheter into the Urethra**

Without Lidocaine anaesthetic  
6ml and 11ml syringes

| Art. No.     | Description   | Volume |
|--------------|---------------|--------|
| <b>CG 6</b>  | Lubricant Gel | 6ml    |
| <b>CG 11</b> | Lubricant Gel | 11ml   |




## Sterile Liquid with 10% Glycerine

**to fill Catheter Balloons**

Thanks to the Glycerine Oil, the balloon size of the catheter remains longer constant

Composition:  
90% Aquadest  
10% Glycerine Oil

| Art. No.        | Glycerine | Filling volume |
|-----------------|-----------|----------------|
| <b>Block 10</b> | 10%       | 10ml           |



## Hydrophilic Urinary Catheter for Paraplegic

**"Do-it-yourself" with water sachet in the pouch**

| Male            |              | Female          |              | Paediatric      |              |
|-----------------|--------------|-----------------|--------------|-----------------|--------------|
| Art. No.        | Size (Fr/CH) | Art. No.        | Size (Fr/CH) | Art. No.        | Size (Fr/CH) |
| <b>ICM 4006</b> | 6            | <b>ICF 2012</b> | 12           | <b>ICP 2006</b> | 6            |
| <b>ICM 4008</b> | 8            | <b>ICF 2014</b> | 14           | <b>ICP 2008</b> | 8            |
| <b>ICM 4010</b> | 10           | <b>ICF 2016</b> | 16           | <b>ICP 2010</b> | 10           |
| <b>ICM 4012</b> | 12           | <b>ICF 2018</b> | 18           |                 |              |
| <b>ICM 4014</b> | 14           |                 |              |                 |              |
| <b>ICM 4016</b> | 16           |                 |              |                 |              |
| <b>ICM 4018</b> | 18           |                 |              |                 |              |



## Hydrophilic Urinary Catheter Kit for Paraplegic

**"Do-it-yourself" KIT – with water sachet in the pouch and attached Urine-Bag**

| Male            |              | Female          |              | Paediatric      |              |
|-----------------|--------------|-----------------|--------------|-----------------|--------------|
| Art. No.        | Size (Fr/CH) | Art. No.        | Size (Fr/CH) | Art. No.        | Size (Fr/CH) |
| <b>IKM 4006</b> | 6            | <b>IKF 2012</b> | 12           | <b>IKP 2006</b> | 6            |
| <b>IKM 4008</b> | 8            | <b>IKF 2014</b> | 14           | <b>IKP 2008</b> | 8            |
| <b>IKM 4010</b> | 10           | <b>IKF 2016</b> | 16           | <b>IKP 2010</b> | 10           |
| <b>IKM 4012</b> | 12           | <b>IKF 2018</b> | 18           |                 |              |
| <b>IKM 4014</b> | 14           |                 |              |                 |              |
| <b>IKM 4016</b> | 16           |                 |              |                 |              |
| <b>IKM 4018</b> | 18           |                 |              |                 |              |

Class IIa device – for single use only

# Intra – Uterine Insemination Catheter

## Artificial fertilization

The injection of sperms into the uterus

Length 17cm



Art. No.  
**IUI**

# Dialysis Catheter


With soft Tip and stepper Dilator  
For venous blood purification



## Straight Dual Lumen

1 Set contents: 1 Flexi-Catheter Dual Lumen,  
1 Stepper Dilator, 2 Injection Caps,  
1 Introducer Needle, 1 Guidewire

| Art. No.            | Size (Fr/CH) | Length (cm) |
|---------------------|--------------|-------------|
| <b>DLCS-K 12-13</b> | 12           | 13          |
| <b>DLCS-K 12-16</b> | 12           | 16          |



## Curved Dual Lumen

1 Set contents: 1 Flexi-Catheter Dual Lumen,  
1 Stepper Dilator, 2 Injection Caps,  
1 Introducer Needle, 1 Guidewire

| Art. No.            | Size (Fr/CH) | Length (cm) |
|---------------------|--------------|-------------|
| <b>DLCC-K 12-13</b> | 12           | 13          |
| <b>DLCC-K 12-16</b> | 12           | 16          |





Marflow AG  
Soodstrasse 57  
8134 Adliswil/Zurich  
Switzerland  
Tel. +41 44 709 01 01  
Fax +41 44 710 21 36  
[marflow@marflow.ch](mailto:marflow@marflow.ch)  
[www.marflow.ch](http://www.marflow.ch)

November 2019

**MarFlow<sup>®</sup>**  
SWITZERLAND