

SafeAir Smoke Evacuator

Effective, silent and easy to use



Effective smoke evacuation

- Instant maximum suction at the tip of the SafeAir Smoke Pencil thanks to powerful vacuum technology
- ULPA filter captures hazardous particles
- Active carbon absorbs smell

Silent

- The vacuum technology gives effective suction with minimal sound

Simple to use

- Automatically synchronised activation with ESU via pop-out sensor
- Intuitive interface, chronological set-up, left to right

SafeAir Smoke Evacuator

Effective, silent and easy to use

Technical specifications

Classification	Class I, type B
Protection class	IP 21, normal tight
Network voltage	220-240VAC - 50 Hz
Power consumption	In operation: 250W Stand-by: 14W
Fuses	2 x T 4AH / 250VAC ceramic fuses
Suction capacity at the tip of the SafeAir Smoke Pencil	High suction: approx. 96l/min Low suction: approx. 45l/min
Vacuum	150 mbar
Level of noise	46 dBA
Filter	ULPA filter with active carbon
Dimensions	450 x 430 x 205 mm
Weight	26.5 kg

Ordering information

EVACUATOR

SFR-0200	SafeAir Smoke Evacuator, incl. one ULPA filter with active carbon	box of 1
----------	---	----------

CONSUMABLES

SFR-FIL-C	ULPA filter with active carbon	box of 1
SFR-FIL	ULPA filter	box of 1
SFR-CON	Tubing to connect ULPA filter and collection bottle	box of 1
SFR-TUB-7	7m extension tubing for SafeAir Smoke Evacuator exhaust	box of 1
SFR-FSW	Footswitch	box of 1
SFR-KIT-1	One-year service kit containing: Air filter, silencer, silicone grease and gloves	box of 1
SFR-KIT-3	3 year service kit containing: 3 lamellas, silencer, O-ring, silicone grease and gloves	box of 1

CE SafeAir Smoke Evacuator

SafeAir AG

D4 Platz 3
CH-6039 Root Längenbold
Switzerland

Email: info@safear.ch
www.safear.ch

VAT: CHE-299.243.601
SafeAir is a trademark of SafeAir
AG Copyright SafeAir 2012

